TAN Bead Nucleic Acid Extractor (Non-sterile)

(Compatible with power 220V-240V only)

Maelstrom 9610 User Manual

039.M1601.X21 Rev V1.0_2022.01

Table of Contents

1.	Introduction	6
2.	Instrument Overview	11
3.	Installation & Get Started	15
4.	User Interface Guidelines	16
5.	Technical Support	38
6.	Cleaning and Maintenance	39
7.	Disposal	39
8.	Patents	40
9.	About Manufacturer	40
10.	Appendix	41

About Manual

The label on the instrument, the User Manual, and other packaging material may contain following symbols:

NEF	
SN	

DEE

Serial number

Catalog number

Manufacturer

The date of manufacture

Keep Away From Sunlight

Non-Sterile

Temp Limit

(Operating)

Humidity Limit

Consult Instructions for Use

The CE mark

In Vitro Diagnostic Medical Device

Authorized European Representative

Caution, consult accompanying documents

Hot surface, risk of burns

Watch your fingers and your hands

The symbol indicates that the product should not be discarded as unsorted waste but must sent to separate collection facilities for recovery and recycling

About Instrument

WARNING

- Use a power cord that meets your country's standard. In case of any questions about the power cord, contact your local distributor or reseller.
- Do not use the instrument with voltage other than the voltage specified on the device.
- Do not use the instrument with damaged power cord or a loose socket.
- When the power plug is disconnected from the AC outlet, please hold the power plug itself instead of pulling the power cord.
- Before working on maintenance, must disconnect the power plug from the outlet in advance.
- Do not pour any liquid on the instrument.
- Do not place any containers with liquid on the instrument. Doing so may cause a fire, electric shock or malfunctions of the instrument.
- Do not touch the power plug or cord if there is a chance of lightning. Failure to observe this may cause electric shocks.
- If you hear thunder or suspect approaching lightning, quickly turn off the power switch and pull the power plug from the AC outlet.
 Failure to observe this may cause a fire or malfunctions.

About Instrument

CAUTION

- Never attempt to remodel the instrument without the permission from the manufacturer. Doing so may lead to a fire or electric shock.
- Do not bump or knock the instrument. Doing so may cause malfunctions.
- Repairs to the instrument must be performed by agencies authorized by Taiwan Advanced Nanotech Inc.
- Only original spare parts from Taiwan Advanced Nanotech Inc can be used on the instrument.
- If the equipment is used in a manner not specified by the manufacturer, the protection given by the instrument may be impaired or invalid.

Maelstrom 9610 is an automated nucleic acid platform designed for high-throughput applications. Specialized spin tips enable super mixing efficiency of magnetic beads and larger processing volume. With intuitive interface and flexible program, Maelstrom 9610 can boost molecular diagnostic laboratory productivity by transforming routine operations into a walk-away solution.

Operation Principle

Maelstrom 9610 uses magnetic rods to collect and transfer magnetic beads and spin tips to mix suspension. Rotating the spin tips will mix suspension very well. Extended magnetic rods can collect magnetic beads efficiently.

Purified nucleic acids can be obtained after cell lysis, nucleic acids absorption, wash and elution.

Intended Use

- Maelstrom 9610 shall be used in combination with TANBead Nucleic Acid Extraction Kit to extract DNA/RNA from specimen.
- Maelstrom 9610 can handle magnetic beads for extraction and purification of nucleic acid.
- This product should be operated by a physician, medical examiner or related professionals.

Environmental Requirements

To avoid shortening of the lifespan of the instrument, use Maelstrom 9610 in a location that meets the following criteria:

- Choose a location that has good air circulation.
- Place Maelstrom 9610 on the table capable of bearing at least the weight of Maelstrom 9610 (≒ 95 kg).
- Do not use Maelstrom 9610 in a location where there is large temperature variability or high humidity.
- Suggested operation condition: Temperature: 10-35°C Relative humidity: 40-80%
- Storage and transportation condition: Temperature: 5-60°C Relative humidity: 20-90%
- Maximum operating altitude: 2000 m

Safety Instructions and Guidelines

- This device can be used with potentially biohazardous materials. Use appropriate personal protective equipment (gloves, safety goggles, lab coat, etc.) for handling and disposing of biohazardous materials.
- Under normal condition, sound pressure level from Maelstrom 9610 does not exceed 80 dB and cause a hazard. Please contact technical support in case of higher sound pressure level.
- This device can be hazardous due to the use of chemical and biohazardous substances.
- Users should adhere to their institutional guidelines for the handling and disposal of all infectious substances used with this device.
- It is important to clean the device after every use. If samples or reagents have been spilled, it is important to clean the instrument immediately to avoid damage or contamination of samples.
- This device is to use with the compatible spin tips. Using incompatible spin tips may cause poor extraction performance.
- Read this user manual in its entirety prior to operating the device. Failure to read, understand, and follow the instructions in the manual may result in damage to the device, injury to laboratory and operating personnel or poor performance.

Safety and EMC Information

Safety Requirements

- The medical device has passed the tests and conformed to the standards of IEC 61010-1:2010, "Safety requirements for electrical equipment for measurement, control, and laboratory use Part 1: General requirements".
- The medical device has passed the tests and conformed to the standards of IEC 61010-2-101:2018, "Safety requirements for electrical equipment for measurement, control and laboratory use
 Part 2-101: Particular requirements for in vitro diagnostic (IVD) medical equipment".

EMC Requirements

- The medical device has passed the tests and conformed to the standards of IEC 61326-1:2012 & IEC 61326-2-6:2012, "Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 2-6: Particular requirements - In vitro diagnostic (IVD) medical equipment".
- The medical device has passed the tests and conformed to following standards:
- EN 61326-1:2013
- EN 61326-2-6:2013

Equipment list

Please note that following accessories accompanying with instrument may vary from region to region.

Power cord x1

Transmission cable x1

Magnetic separator x1

Tips combo pack x1

Auto Tube base x 6

=Remark=

Power cord varies from region to region.

Auto Tube base is required in use of Auto Tube reagent kits.

User manual x1

ITEM	SPECIFICATION
Model	Maelstrom 9610
Weight	Approx. 95 kg
Dimensions	870(W) x 575(L) x 700(H) mm
Power rating	220-240Vac, 50/60Hz, 3.5A (Class I)
Fuse	250V, 5A
Max. Throughput	96 samples per run
Process. volume	50µl -1,600µl
Spin speed	up to 3000 rpm
Heater	4 independent heating plates
Magnetic rod	>3900 gauss
Display	7-inch touchscreen
UV & HEPA	UV: UV-C type 4W, HEPA: E10 Class

Please note that this machine is compatible with 220-240V power only. Using wrong power may lead to malfunction or damage to the instrument.

Rear View

- 1 USB type A for use of upgrading HMI program
- 2 USB type A for use of upgrading PLC program
- 3 Power inlet
- $\textcircled{4} \mathsf{Fuse}$

Please note that only TANBead technicians and authorized technicians or service engineers shall engage in HMI and PLC program upgrade.

Control Panel

	ia
Taiwab Nanotech Nanotech	Buzzer
Tip Run Report	
UV/HEPA Program Setting The Setting	Pause (
	(💿 Stop

To turn off the alarm

To pause the program. Press again to resume the program.

To turn on/off the LED light.

To initiate the program

Emergency to stop machine

3. Installation & Get Started

Installing Maelstrom 9610 is easy. Please follow the guide below to correctly set up the instrument.

Please note that this instrument weights around 95kg, it is highly recommended to have 3 or more persons to handle when taking it out from the box and must pay attention to the safety.

Step 1:

Take out the instrument from the box (3 or more persons to handle is highly advised)

Step 2: Remove the cushions around the instrument

Step 3:

Place the instrument to a flat-table that can bear 120kg weight

Step 4:

Connect the power. Please note that this instrument is compatible with 220-240V power only. Using wrong power source will lead to malfunction or damage to the instrument

Step 5:

Power on the instrument and it will perform initialization. Once initialization is done, please remove Y-axis protector from the instrument and the instrument is ready to use.

Y-axis protector

Create an User ID

An user ID is required to access the instrument. For a brand new machine, there is a need to create an user in advance. Click the icon to create your own user ID.

СК	197	5									
ESC				4		6					•
CLR		w	E		т		U		C		
									К		
								N.		Ente	

 Use the keyboard to input user ID and then click [Enter] to save and exit

User ID: max. 12 characters

Save User ID

After creating the user ID, please must click \bigcirc to save and then program will automatically lead you access the system.

Switch user

To switch user, there is no need to power off/on the instrument. Take advantage of icon to switch user quickly.

Possible Error Message

Some messages may occur when you try to log in or create an user ID.

Main Page

There are six functions that can be used on main page.

To test on mounting and ejecting spin tips

To execute the extraction program

To export extraction records

To use UV & HEPA to clean the instrument

To edit, add, and delete extraction program

To set basic parameter of instrument

Tip (Test on Mount and Eject Spin Tips)

- □ Tip Point (Setting): Parameter setting value upon mounting tips
- □ Tip Point (Actual): Actual tip-mounting parameter value
- **D** Eject Tip Point (Setting): Parameter setting value upon ejecting tips
- □ Eject Tip Point (Actual): Actual tip-ejection parameter

=Remark= Actual value needs to equal to Setting value

How to test mounting and ejecting spin tips

Step 1: Click "All Home" to set X/Y/Z-axis to zero

Step 2: Place tip combo pack (for use of mounting spin tips) or blank 96DW plate (for use of ejecting spin tips) to #8 plate position

Tip function is designed to test the performance of mounting and ejecting spin tips.

Step 3: Use "▶" and "◄" button to move #8 plate underneath the extraction module

Step 4: Click "Mount" icon to execute mounting spin tips (tip combo pack required)or click "Eject" icon to execute ejecting spin tips (blank 96DW plate required)

Step 5: Check if all spin tips are mounted to spin shaft without evident gap or if all spin tips are completely ejected from spin shaft

Run (Execution of Extraction Program)

Run function is designed to execute nucleic acid extraction program stored in the instrument.

Maelstrom 9610 is capable of storing up to 25 programs. Please follow the guideline of reagent kits and correctly select the program.

Step 1: Click the program first, and then click [OK] to continue

Step 2: Click [Next] to continue

Run (Execution of Extraction Program)

Step 3: Click space to input "Kit Lot" and "Tip Lot", if need. Click [Next] to continue. Both Kit Lot and Tip Lot can bear 18 characters input.

Step 4: Input sample ID, if need. Press the circle first, and then use onscreen keyboard or USB barcode scanner to scan in if 1-D barcode is available. Sample ID can bear 18 characters input. Once completed, Press [Next] to continue.

Pos.: Hole position on the plate
No.: Total number of Sample ID
ID: Name of Sample ID

=Remark=

Maelstrom 9610 can be compatible only with USB-type infrared barcode scanner.

Run (Execution of Extraction Program)

Step 5: Follow the guide on screen to place the reagent plate onto correct plate accordingly. Please note that the priority well plate must be spin tips compo pack, please place the eighth plate.

=Remark=

Eject tip plate is allowed to selected from NO.2 to NO.5 plate according to user's needs.

Step 6: Ready to run the extraction program. Click [Yes] to begin the operation. If you choose [No], the process will go back to previous page.

During the Run

Understanding every symbol in the operation process is important. It helps you clearly know the status of nucleic acid extraction.

- Start time: Time of extraction program started
- Est. End time: Estimated time the extraction program coming to the end
- □ Speed: Spin speed
- □ Volume: Liquid volume in the working well

During the Run

Understanding every symbol in the operation process is important. It helps you clearly know the status of nucleic acid extraction.

Mixing

Beads collection

In the course of beads collection, whether you click Pause button or Stop button, the instrument will not act immediately until the collection process comes to the end.

Vaporing EtOH

Before entering Elute process, the protocol will must work on Vapor process to evaporate EtOH. Poor evaporation on EtOH will influence the qPCR data significantly.

Program finished

When a run comes to the end, the Buzzer will ring up few seconds.

Resume a run

During the operation, if you click Pause to temporarily stop the program, you can click Resume to continue the operation.

Report (Export extraction records)

Report function enables you to view experiment records and the identity that performed the experiment.

▲ Maelstrom 9610 is designed to store the last 15 experiment records.

ľ			2019/02/16	00:10
0	Риссист	Position	Sample ID	
	F E U Y E AM	A01	ABCDEFGHIJKLMNOPQR	
	ABCDEFGHIJ	A02	ABCDEFGHIJKLMNOPQR	
	Kit Lat	AOB	ABCDEFGHIJKLMNOPQR	
	NIL LUL	A 0 4	ABCDEFGHIJKLMNOPQR	
	ABCDEFGHIJKLMNOPOR	A05	ABCDEFGHIJKLMNOPOR	
		A06	ABCDEFGHIJKLMNOPQR	
	Tip Lot	A07	ABCDEFGHIJKLMNOPQR	
	ABODEEGHTIKIMNOPOP	A 0 8	ABCDEFGHIJKLMNOPQR	
\leftarrow	Abeber offisikennor on	A09	ABCDEFGHIJKLMNOPQR	
		AÍO	ABCDEFGHIJKLMNOPQR	
		A11	ABCDEFGHIJKLMNOPQR	-
		A12	ABCDEFGHIJKLMNOPQR	

By clicking View button, you can navigate the detail of sample ID.

UV / HEPA

It is highly advised to use UV to clean the instrument after each experiment.

Click HEPA icon to switch to HEPA function Same with UV, set HEPA time first and then click icon to start

Before turn on UV, please must close the cover door of turntable. The recommended lighting time is at least 30min. The more, the better.

Program (Edit/Delete/Add a Program)

Program function enables you to customize the extraction protocol.

Click a protocol first, and then choose a function to proceed

ľ	Edit can be used to adjust the setting of existing protocols, which contains well volume, preheat temp, spin direction, extraction steps, time of mixing and collect.
(<u>)</u>	To add a new program, you need to use Copy to duplicate a program first, and then use Edit to write the protocol.
×	Delete a program no longer used to save the memory space. Maelstrom 9610 can store up to 25 programs and each program can support totally 35 steps.

Edit a Program

Ľ							2021/0	2/25 11	: 08
	Ш	1	₽[ABCDEF	HIJ				20
	Plate	1	2	3	4	5	6	7	8
	Volume (ul)	1600	1600	1600	1600	1600	1600	1600	
	Preheat	100	100	100			100	\sim	
	Action).).		For./Rev.	For.	For./Rev. U/D	For.U/D		
	Name	ABCD	ABCD	ABCD	ABCD	ABCD	ABCD	ABCD	ABC
				N	ext				

A Press existing value to change the setting

Program	Name of program, up to 10 characters
Plate	No. of plate
Vol.(µl)	Total processing volume
Preheat (°C)	Heating a plate to a particular temp before mixing step Preheat temp range from 25°C to 100°C
Action	Way of binding during the mixing process For.: Forward For./Rev.: Forward and reverse (spin direction will change regularly) For. U/D: Forward with up-and-down move For./Rev. U/D: Forward and reverse with up-and-down move
Name	Name of plate, which following buffer may be involved in LB: Lysis buffer WB1: Wash buffer 1 WB2: Wash buffer 2 MB: Magnetic beads EB: Elute buffer

Please note that name of plate 8 is unalterable.

Edit a Program

ľ	📰 10	EFGHIJ	<u>[₁</u>] ABC	[2] ABC	<u>[3]</u> [4] ABC ABC	2019/02 [5] [6] [ABC ABC AB	/16 00:10 7] <u>[8]</u> 8C TIP
	Step Plate	e Temp.	Mixing (min)	Spin (rpm)	Collect (sec)	Vapor Pause (min)	ŝ.,
	1 1	130	999.9	3000	999	999.9 O off)
	2 2	130	999.9	3000	999	999.9 OFF)
5	3 3	130	999.9	3000	999	999.9 Ooff)
_	4 4		999.9	3000	999	999.9	
	5 5		999.9	3000	999	999.9	•

Step	Extraction step Maelstrom 9610 supports up to 35 steps in a program
Plate	No. of plate
Temp	Please adjust to appropriate protocol temperature before running the test which is suitable for different reagents. Protocol temp and Preheat temp are different. When executing a program, preheat temp will be fulfilled in advance, and protocol temp will be acted as the target plate is about to be used.
Mixing	Time of mixing, range from 0 to 999.9 mins
Spin	Speed of spin (mixing), range from 500 to 3000 rpm
Collect	Time of collecting beads, range from 0 to 999.9 sec
Vapor	Time of evaporate EtOH, range from 0 to 999.9 mins
Pause	Will automatically pause after executing the current step. If additional reagents need to be added during the extraction process, this step can be taken.

Please note that name of plate 1 and plate 8 is unalterable.

Add a new program

- 1. Touch one of standard programs button.
- 2. Press Copy button and fill in blank space.
- 3. Choose new program and press Edit button to set the relevant parameters.

▲ Touch one of standard programs button and press Copy button.

Confirm the number to add the program. The keyboard will appear after press on destination number. Once complete, click Yes to add the program.

Delete a Program

To delete a program, press the program no longer need and then click Delete icon to execute.

▲ Press the unwanted program, and press the Delete icon

Click Yes to confirm the deletion

Setting

Setting is designed for user to change basic system parameter setting, such as display language, time info and barcode scanner management.

Admin user is allow to assess Language, Time, and Barcode functions.

Language

Maelstrom 9610 supports multi-language display, including English, Traditional Chinese, Simple Chinese, Korean and Spanish.

Time

When a new instrument is delivered, please adjust to the right time zone.

Barcode

Scan code to input related information. e.g. lot number, sample ID.

Setting (Language)

Maelstrom 9610 supports multi-language operating interface including English, Traditional Chinese, Simple Chinese, Korean and Spanish.

▲ Press Language icon to access

مر Settings	2021/07/08 17:33:0
Language	
	English 繁體中文 简体中文
Previous Page	Español 한국어
Main Page	

▲ Choose the different languages

Setting (Time)

Set correct time for the system according to the local time

Click icon to set Date and Time for the system; once set up is done, press [OK] to apply the latest setting.
4. User Interface Guidelines

Setting (Barcode)

Maelstrom 9610 can be compatible with USB-type infrared barcode scanner.

▲ Click the icon to ON/OFF of barcode scnner

If barcode scanner is set OFF, even USB-type barcode scanner is connected to Maelstrom 9610, scanned data will not be presented.

5. Technical Support

Taiwan Advance Nanotech Inc. provide post-sales service and technical support. In case of any questions, please try to contact our authorized distributor nearest to you or our company for more help.

If there is a need to contact us, you can call our rep number at +886-3-3167568 or email us at service@tanbead.com for further assistance.

To efficiently solve the problems and answer your questions, please provide us your instrument serial number when you talk to our technician personnel. Thanks for your cooperation.

6. Cleaning and Maintenance

- It is important to clean the device after every use. If samples or reagents have been spilled, it is important to clean the device immediately to avoid damage or contamination of samples.
- Wear gloves and appropriate personal protective equipment. If the device is used with biohazardous materials, dispose of any cleaning materials used in accordance with your institutional guidelines.
- The device may go through a run with the magnetic rods unprotected. If this happens, the magnetic rod needs to be cleaned immediately.
- To clean the magnetic rods, wipe with a soft cloth dampened with pure water. Do not use alcohol solvent.
- If the magnetic rods cannot be cleaned, please contact TANBead (Tel: +886-3-3167568) for technical assistance.

7. Disposal

In case the product is to be disposed of, the relevant legal regulations are to be observed. As disposal regulations may differ from country to country, please contact your supplier if necessary. Please follow your institutional and country-specific requirements to handle the disposal of accessories. Device must be decontaminated prior to disposal.

ENGLISH

8. Patent

Patent List	
USA	US09616398B2
EU	EP2937136
Canada	CA2862946
Japan	JP6151735B2
Korea	KR101696517B1
China	CN104971638B
Taiwan	TWI526245B
WIPO	WO2016127292

9. About Manufacturer

Here comes the detail of instrument manufacturer and its EC Rep.

- Manufacturer : Taiwan Advanced Nanotech Inc.
- Manufacturer Add : 4F, No. 188, Wenhe Rd., Guishan Dist., Taoyuan City 333, Taiwan (R.O.C.)
- Manufacturer Tel: +886-3-3167568

EC REP

- Contact details : mdi Europa GmbH
- Address : Langenhagener Street 71, D-30855 Langenhagen,Germany
- Telephone: +49-511-3908-9530

Fuse Replacement on Maelstrom 9610

Fuse Specification :

- UL recognized component
- Rated 5A, 250VAC x2pcs
- UFE, Size 5.20 x 20mm

Safety precautions of replacing fuses

The following safety precautions will prevent injury to personnel and damage to equipment. These are the MINIMUM safety precautions for replacing fuse.

- 1. Be sure the power is off in the circuit and the circuit is discharged before replacing a fuse.
- 2. Use an identical replacement fuse.
- 3. Remove any corrosion from the fuseholder before replacing the fuse.
- 4. Be certain the fuse properly fits the fuseholder.

How to fix the fuse

- 1. Turn the power off and discharge the circuit before removing a blown fuse.
- 2. Use a fusepuller (an insulated tool) when you remove a fuse from a clip-type fuseholder.
- 3. Remember, always replace the faulty fuse with one equal in size, amperage, rating and type.
- 4. Turn the main power switch back on and see if electricity is restored.

How to use Maelstrom Master to update extraction protocols in the instrument?

Connect the instrument to your PC

Power on the instrument and then use the cable (USB A to B) to link the instrument to your PC (Type A to PC, Type B to instrument).

Execute the utility program "Maelstrom 9610 Master"

You can find the utility program from the CD that comes with the instrument. Load the CD to the PC and use Windows Explorer to find the folder named "Maelstrom Master". Then, double click the utility "TANBead Maelstrom 9610 Master".

The utility program will automatically detect the connection between the instrument and PC. If successfully connected, the pop-out window will tell.

How to use Maelstrom Master to update extraction protocols in the instrument?

Program Editing

After connection, the utility program will automatically read all protocols stored in the instrument.

Select the protocol first and choose the function to performed. Please note that the protocol you select will be marked in Blue.

How to use Maelstrom Master to update extraction protocols in the instrument?

Remember to SAVE settings

After adjustment of protocol contents, please MUST click "Save" to get most updated contents effective. Once you click "Save", the utility program will automatically update the protocol contents to the instrument.

🖳 TAN Bea	d Maelstrom 9610	Maste	r						_	.	
Program	Ē	24	ļ						± S	ave	4
T		63	85-N	/lodifi						-1	\mathbf{X}
Edit		Pla	te	Temp	Mixing (min)	Spin (rpm)	Collect (sec)	Vapor (min)	Paus	se	
	1	3	•	30	1	3000	30	0	OFF	•	^
	2	2	•	55	1	3000	0	0	OFF	•	
	3	1	•	55	20	3000	0	0	OFF	•	
\leftarrow	4	2	-	55	0	3000	30	0	OFF	•	
Previous	5	1	•	55	10	3000	30	0	OFF	•	
Page	6	2	-	55	2	3000	30	0	OFF	•	
	7	3	•	30	2	3000	30	0	OFF	•	
	8	4	•	0	2	3000	30	0	OFF	٠	~

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