



Taiwan  
Advanced  
Nanotech



TAN Bead Nucleic Acid Extractor  
(Non-Sterile)

039.L2601.X01




















# Maelstrom 9610 LH User Manual

# Contents

1.	Introduction .....	6
2.	Instrument Overview .....	12
3.	Installation & Get Started .....	14
4.	Command List .....	15
5.	Parse "GetErrorCode" Response .....	31
6.	Response Error Code .....	32
7.	Technical Support .....	33
8.	Cleaning and Maintenance .....	34
9.	Disposal .....	34
10.	Patent .....	35
11.	About Manufacturer .....	35
12.	Advisory Notices .....	35
13.	TANBead Instrument Warranty Policy.....	36

# About Manual

The label on the instrument, the User Manual, and other packaging material may contain following symbols:

- |   |  |  |                                    |
|---|--|--|------------------------------------|
|                    | Catalog number   |   | Temp Limit (Operating)             |
|                    | Serial number  |   | Keep Dry                           |
|                    | Manufacturer   |   | Keep Away From Sunlight            |
|                    | The date of manufacture  |   | EU authorized representative       |
|                    | The CE mark  |   | In Vitro Diagnostic Medical Device |
| <br>eIFU Indicator | Instructions for Use   |   | Humidity limit                     |
|                    | Non-Sterile  |   | Caution                            |
|                  | Biologic Risks   |  | Biologic risks                     |
|                  | Hot surface, risk of burns   |  |                                    |
|                  | Watch your fingers and hands   |  |                                    |
|                  | The WEEE symbol, indicating separate collection for WEEE- Waste of electrical and electronic equipment |  |                                    |

## About Instrument

### Warning

- Use a power cord that meets your country's standard. In case of any questions, contact your local distributor for assistance.
- Maelstrom 9610 LH operates within the voltage range of 100 Volts and 240 Volts.
- Do not use the instrument with damaged power cord or a loose socket.
- To disconnect the power plug from the AC outlet, hold the power plug itself instead of pulling the power cord.
- Prior to performing the maintenance, make sure to disconnect the power plug from the outlet.
- Do not pour any liquid on the instrument.
- Do not place any containers with liquid on the instrument. Doing so may cause a fire, an electric shock or malfunctions of the instrument.
- Do not touch the power plug or cord if there is a chance of lightning. Failure to observe this may cause electric shocks.
- If you hear a thunder or suspect an approaching lightning when in use, turn off the power switch and disconnect the power plug from the AC outlet immediately. Failure to observe this may cause a fire or malfunctions.

# About Instrument

## Caution

- Never attempt to remodel the instrument without the permission from the manufacturer. Doing so may lead to a fire or an electric shock.
- Do not subject the instruments to any impacts and do not knock it. Doing so may cause malfunctions.
- Any repairs to the instrument must be performed by agencies authorized by Taiwan Advanced Nanotech Inc.
- Only use the original spare parts supplied by Taiwan Advanced Nanotech Inc on the instrument.
- If the equipment is used in a manner not specified by the manufacturer, the protection given by the instrument may be impaired or invalid.

# 1. Introduction

## About

Maelstrom 9610 LH is easily integrated into Liquid Handling Workstation, introducing a new tip sensor which automatically reminds user that the instrument is failure picking up or leaving tips.

## Principle

Maelstrom 9610 LH uses patented magnetic beads, spin tips and reagent kits technologies , consisting of three major process: isolation, purification and concentration.

## Intended Purpose/ Intended Use

The Maelstrom 9610 LH is intended for medical laboratory use by trained personnel in automation environments. The instrument is intended for automated transfer and processing of magnetic particles in a microplate format, to extract and purified nucleic acids from human samples. The purified nucleic acid can be used with any downstream application employing PCR-based qualitative, semi-quantitative and quantitative assays.

# 1. Introduction

## Environmental Requirements

To avoid shortening the lifespan of the instrument, use Maelstrom 9610 LH in a location that meets the following criteria:

- Choose a location with good air circulation.
- Place Maelstrom 9610 LH on the table that can bear at least 30 kg
- Do not use Maelstrom 9610 LH in a location where is with huge temperature and humidity variability.
- Operate condition:  
Temperature: 10-40°C  
Relative humidity: 40-80%
- Storage and transport condition:  
Temperature: 5-50°C  
Relative humidity: 20-85%
- Maximum operate altitude:  
2000m

# 1. Introduction

## Safety Instructions and Guidelines

- This device can be used with potentially biohazardous materials. Use appropriate personal protective equipment (gloves, safety goggles, lab coat, etc.) for handling and disposing of biohazardous materials.
- Under a normal condition, sound pressure level from Maelstrom 9610LH does not exceed 80dB and does not cause a hazard. Please contact technical support for assistance in case of a higher sound pressure level.
- This device can be hazardous due to the use of chemical and biohazardous substances.
- Users should adhere to their institutional guidelines for the handling and disposal of all infectious substances used with this device.
- It is important to clean the device after every use. If samples or reagents have been spilled, clean the instrument immediately to avoid damage or contamination of samples.
- This device is to use with the compatible spin tips. Using incompatible spin tips may cause poor extraction performance.
- Read this user manual completely prior to operating the device. Failure to read, understand, and follow the instructions in the manual may result in damage to the device, injury to laboratory and operating personnel or poor performance.
- If any serious incident occurs, please report to the manufacturer and the competent authority of the member state in which the user and/or the patient is established.

# 1. Introduction

## Safety Requirements

- The device has passed the tests and conformed to the standards of IEC 61010-1:2010+A1:2016 (Edition 3.1) and/ or EN 60101:2010+A1:2019, "Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 1: General requirements".
- The device has passed the tests and conformed to the standards of IEC 61010-2-101:2018 with IEC 61010-1:2010 + A1:2016 and/ or EN 61010-2-101:2017 with EN 61010-1:2010 + A1:2019, "Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-101: Particular requirements for in vitro diagnostic (IVD) medical equipment".

## EMC Requirements

- The device has passed the tests and conformed to the standards of IEC 61326-1:2020 / EN IEC 61326-1:2021 & IEC 61326-2-6:2020 / EN IEC 61326-2-6:2021, "Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 2-6: Particular requirements - In vitro diagnostic (IVD) medical equipment".

# 1. Introduction

## Accessory


Following accessories may vary region-to-region.



Power cord



Power adapter




Magnetic Base



Transmission cable

## Consumables

Maelstrom 9610LH uses specially designed consumables for optimal processing. Use of other types of plates may damage the instrument and compromise the warranty


Item	Small Magnetic Rod (B)
<b>96 Spin Tips Assembled Box</b>	(Ø2.2) Assembled Box 
<b>96 DWP</b>	96 Deep Well Plate 
<b>Auto tube</b>	16 auto tube  16 base B 

## 2. Instrument Overview



ITEM	SPECIFICATION
Model	Maelstrom 9610 LH
	Small Magnetic Rod (B)
Magnetic Rod Diameter	Ø2.2 mm
Weight	Approx.25kg
Dimensions	542(W) x 285(D) x 305(H)mm
Power rating	I/P: 100~240V AC, 3A Max, 50/60Hz O/P: DC 24V, 8.3 A
Max. Throughput	6 to 96 samples/run
Process volume	50~1,600µL
Spin speed	500~3,000rpm
Heater	1 independent heating plate
Heating Plate	U bottom
Temperate Range	RT to 100°C
Magnetic rod	>3,900 gauss

## 2. Instrument Overview



- ① USB type A for connecting with PC
- ② LED indicating Power ON/OFF
- ③ Power switch
- ④ Power inlet


### Connections :

**Computer minimum requirements:**

**Processor:** Intel Core i5, 2.50GHz

**RAM:** 8 GB

**OS:** Windows 10



## 3. Installation & Get Started

Please note that this instrument weight is around 30kg, it is highly recommended to have 2 individuals or above to handle when taking it out from the box and be sure to work safety.

### Step 1:

Take out the instrument from the box, 2 or more individuals is highly recommended.

### Step 2:

Remove the cushions around the instrument.

### Step 3:

Place the instrument on a flat table that bears over 30 kilograms.

### Step 4:

Connect the power. Please note that this instrument is compatible with AC 100-240V power only. Using wrong power source will lead to malfunction or damage.

### Step 5:

Connect the system to the computer (Window 10 OS, RAM >8G, Intel core i5 or above) through USB port. Install the driver on the computer, which is provided by the manufacturer.

### Step 6:

Power on the instrument and it will perform initialization.

### Step 7:

Operation protocols may vary between reagents. Please follow the protocols listed in the introduction manual for each reagent kit. General procedures will include:

- a. Place the "Spin Tips Assembled Box" onto the heating block.
- b. Mount the "Spin Tips" and replace it with the "Deep Well Plate" containing corresponding reagent.
- c. Perform the following extraction procedure using the magnetic beads in the "Deep Well Plate" for transfer.
- d. Replace the "Deep Well Plate" with other plates as needed to proceed.

## 4. Command List

### Link Control Protocol

USB COM Port = COM50 (com0com-serial port emulators)

Baud Rate = 115200

Data Bit = 8

Stop Bit = 1

Parity Bit = None

By using \n – It prints a new line

### WELL

	Plate Position	Volume	Action	Placeholder	Placeholder
WELL	2	{Volume}	{Action}	0	0

WELL function : Set up the volume and action

{Volume} : Specified as microliter range from 50 to 1600

{Action Type} : {0} Forward

{1} Reverse

{2} Forward U/D

{3} Reverse U/D

For example :

Send to M9610LH :

"WELL 2 800 0 0 0\n" {Volume=800 $\mu$ L} · {Action=0}

M9610LH Response

"WELL ACK\n" To respond after command received

"WELL OK\n" To respond after command accomplished

Any optional field for which there is no meaningful data contains a zero (0) or a space as a placeholder.

## 4. Command List

### PMOV

	Plate Position
PMOV	{Plate Position}

PMOV function : Set up the plate position

{Plate Position} : {1} Move to left position (only using for 48 series)  
{2} Move to right position  
{3} Move to exchange position  
{4} Distance from initial point 2 mm

For example :

Send to M9610LH :

"PMOV 2\n"                      {Plate Position = 2}

M9610LH Response :

"PMOV ACK\n"

To respond after command received.

"GO\_FRONT OK\r\n"

To respond after shield plate moving to front limit position.

"PMOV OK\n"

To respond after command accomplished.

## 4. Command List

### PMIX

	Plate Position	Spin Time	Spin Speed
PMIX	2	{Spin Time}	{Spin Speed}

PMIX function : Set up spin time and speed

{Spin Time} : Specified as minutes from 0.1 to 60

{Spin Speed} : Specified as RPM from 500 to 3,000

For example :

Send to M9610LH :

"PMIX 2 0.5 3000\n"

{Spin Time = 0.5(min)}

{Spin Speed = 3000(rpm)}

M9610LH Response :

"PMIX ACK\n"

To respond after command received.

"GO\_REAR OK\r\n"

To respond after shield plate moving to rear limit position.

"PMIX GO\n"

To respond after moving to the position and start the real mix action.

"GO\_FRONT OK\r\n"

The mix time would start to count down.

To respond after shield plate moving to front limit position.

"GO\_REAR OK\r\n"

To respond after shield plate moving to rear limit position.

" PMIX OK\n"

To respond after command accomplished.

## 4. Command List

### PCOL

	Plate Position	Beads Collection
PCOL	2	{Collect Time}

PCOL function : Set up the collective time

{Collect Time} : Specified as minutes from 0.1 to 60

For example :

Send to M9610LH :

"PCOL 2 0.5\n"                    {Collect time = 0.5(min)}

M9610LH Response :

"PCOL ACK\n"

To respond after command received.

"GO\_REAR OK\r\n"

To respond after shield plate moving to rear limit position

"PCOL POS1\n"

To respond after moving to the position 1 and start the real collect action.

"PCOL POS2\n"

The collect time would start to count down. To respond after moving to the position 2 and start the real collect action.

"PCOL POS3\n"

The collect time would start to count down. To respond after moving to the position 3 and start the real collect action.

"PCOL POS4\n"

The collect time would start to count down. To respond after moving to the position 4 and start the real collect action.

"PCOL POS5\n"

The collect time would start to count down. To respond after moving to the position 5 and start the real collect action.

"GO\_FRONT OK\r\n"

The collect time would start to count down. To respond after shield plate moving to front limit position.

"GO\_REAR OK\r\n"

To respond after shield plate moving to rear limit position.

"PCOL OK\n"

To respond after command accomplished.

## 4. Command List

### PVAPOR

	Plate Position	Vapor time
PVAPOR	2	{Vapor Time}

PVAPOR Function : Set up the vapor time

{Vapor Time} : Specified as minutes from 0.1 to 60

For example :

Send to M9610LH :

"PVAPOR 2 0.5\r\n"

{Vapor time=0.5(min)}

M9610LH Response:

"PVAPOR ACK\r\n"

To respond after command received.

"GO\_FRONT OK\r\n"

To respond after shield plate moving to front limit position and start the real vapor action. The vapor time would start to count down.

"GO\_REAR OK\r\n"

To respond after shield plate moving to rear limit position.

"GO\_FRONT OK\r\n"

To respond after shield plate moving to front limit position.

"PVAPOR OK\r\n"

To respond after command accomplished.

## 4. Command List

### TMNT

	Plate Position
TMNT	2

TMNT function : Pick up tips  
For example:

Send to M9610LH :  
"TMNT 2\n"

M9610LH Response :

"TMNT ACK\n"

"GO\_REAR OK\r\n"

"GO\_FRONT OK\r\n"

"TMNT OK\n"

To respond after command received.

To respond after shield plate moving to rear limit position.

To respond after shield plate moving to front limit position.

To respond after command accomplished.

## 4. Command List

### TEJT

	Plate Position
TEJT	2

TEJT function : Leave tips

For example :

Send to M9610LH :

"TEJT 2\n"

M9610LH Response :

"TEJT ACK\n"

"GO\_REAR OK\r\n"

"GO\_FRONT OK\r\n"

"TEJT OK\n"

To respond after command received.

To respond after shield plate moving to rear limit position.

To respond after shield plate moving to front limit position.

To respond after command accomplished.

## 4. Command List

### TMPE

	Placeholder	Temperature
TMPE	0	{Temperature}

TMPE function : Set up the temperature

{Temperature} : Specified as degrees Celsius from RT to 130

For example :

Send to M9610LH :

"TMPE 0 40\n" {Temperature=40 (deg C)}

M9610LH Response :

"TMPE ACK\n" To respond after receiving command and start the heating process if necessary.

"TMPE OK\n" To respond after temperature reaches the target or after turning off the heating process.

Caution: Please turn off TMPE after heating process.

Send to M9610LH :

"TMPE 0 0\n" {Temperature = 0 (turn Off)}

M9610LH Response :

"TMPE ACK\n" To respond after receiving command and start the heating process if necessary.

"TMPE OK\n" To respond after temperature reaches the target or after turning off the heating process.

## 4. Command List

### PAUS

	Pause or Resume
PAUS	{Pause or Resume}

PAUSE function : PAUSE is only using during spinning period

{Pause or Resume} : {0} Resume  
 {1} Pause

For example :

Send to M9610LH :  
 "PMIX 2 0.5 3000\n"

M9610LH Response :

"PMIX ACK\n"

To respond after command received.

"GO\_REAR OK\r\n"

To respond after shield plate moving to rear limit position

"PMIX GO\n"

To respond after moving to the position and start the real mix action. The mix time would start to count down.

Send to M9610LH :

"PAUS 1\n"

pause mix command

M9610LH Response :

"PAUS ACK\n"

To respond after command received. The mix time would pause to count down.

"PAUS OK\n"

To respond after command accomplished.

"GO\_FRONT OK\r\n"

To respond after shield plate moving to front limit position

## 4. Command List

### PAUS

Send to M9610LH :  
"PAUS 0\n"

resume mix command

M9610LH Response :  
"PAUS ACK\n"  
"PAUS OK\n"  
"GO\_REAR OK\n"

To respond after command received.  
To respond after command accomplished.  
To respond after shield plate moving to rear limit position.

"PMIX GO\n"

To respond after moving to the position and start the real mix action. The mix time would resume to count down.

"PMIX OK\n"

To respond after command accomplished.

### STOP

STOP function: The instrument will be reset and disconnected. And this function is only used for spinning period and idle.

For example :  
Send to M9610LH :  
"STOP\n"

M9610LH Response :  
"STOP ACK\n"  
"STOP OK\n"

To respond after command received.  
To respond after command accomplished.

p.s. When the instrument is executing a STOP order, the stop process needs to be accomplished. Otherwise, the operator can't give a command.

### ORIGIN

ORIGIN function: Three axis (X, Y and Z) moving to basis point

For example :  
Send to M9610LH :  
"ORIGIN\n"

M9610LH Response :  
"ORIGIN ACK\n"  
"ORIGIN OK\n"

To respond after command received.  
To respond after command accomplished.

## 4. Command List

### Read Temp

ReadTemp function : Reading temperature

For example :

Send to M9610LH :  
"ReadTemp\n"

M9610LH Response :  
"ReadTemp 0.00 19.79\n"

To respond after command received.

"ReadTemp {Placeholder (48LH only)} {Temperature}"

### GetFWVer

GetFWVer function: Getting firmware version

For example :

Send to M9610LH :  
"GetFWVer\n"

M9610LH Response :  
"GetFWVer V1.0.0.2T1\_9610LH\n"

To respond after command received.

"GetFWVer {FW Version}"



## 4. Command List

### TMPR

Send to M9610LH :  
"ReadTemp\n"                      Read temperature command

M9610LH Response :  
"ReadTemp 478.55 44.40\n"  
To respond after command received.  
{Temperature= 44.40}

Send to M9610LH :  
"TMPR 0\n"                      Close TMPR

M9610LH Response :  
"TMPR ACK\n"                      To respond after command received.  
"TMPR OK\n"                      To respond after command accomplished.

### GetErrorCode

GetErrorCode function: Getting error code

For example :  
Send to M9610LH :  
"GetErrorCode\n"

M9610LH Response :  
"GetErrorCode 1111111 00000030\n"  
To respond after command received.  
"GetErrorCode {state code(7 ASCII)} {error code(8 ASCII)}"

## 5. Parse "GetErrorCode" Response

1<sup>st</sup> argument: State code (7 ASCII)

1	2	3	4	5	6	7
---	---	---	---	---	---	---

### Func 1 Plate state

- 1) Basis point
- 2) Remove position
- 3) Left position
- 4) Right position
- 5) Moving
- 6) Abnormal

### Func 2 Mixing state

- 1) Normal
- 2) Spinning
- 3) Pause
- 4) Abnormal

### Func 3 Collection state

- 1) Normal
- 2) Elution
- 3) Collection
- 4) Abnormal

### Func 4 Heater state A (position2)

- 1) Stop heating
- 2) Reach setting temperature
- 3) Heating up
- 4) Placeholder
- 5) Abnormal

### Func 5 Heater state B (position1) Func 6 Tip state

- 1) Stop heating
- 2) Reach setting temperature
- 3) Heating up
- 4) Placeholder
- 5) Abnormal

- 1) Normal
- 2) Mounting
- 3) Ejecting
- 4) Abnormal

## 5. Parse "GetErrorCode" Response

1<sup>st</sup> argument: State code (7 ASCII)

1	2	3	4	5	6	7
---	---	---	---	---	---	---

**Func 7** Vapor state

- 1) Normal
- 2) Vapor
- 3) Abnormal

## 5. Parse "GetErrorCode" Response

2<sup>nd</sup> argument: Error code (8 ASCII)

A	B	C	D	E	F	G	H
---	---	---	---	---	---	---	---

**Func A** X-axis motor

- 0) Normal
- 1) Placeholder
- 2) Limit deviate
- 3) Calibration failed

**Func B** Y-axis motor

- 0) Normal
- 1) Placeholder
- 2) Placeholder
- 3) Calibration failed

**Func C** Z-axis motor

- 0) Normal
- 1) Placeholder
- 2) Placeholder
- 3) Calibration failed

**Func D** Slew drive motor

- 0) Normal
- 1) Unusual rotational speed

**Func E** Drip-proof motor

- 0) Normal
- 1) Unusual front position
- 2) Unusual rear position

**Func F** Heater A (position2)

- 0) Normal
- 1) Over heating
- 2) Unusual heating time
- 3) Temperature issue

## 5. Parse "GetErrorCode" Response

2<sup>nd</sup> argument: Error code (8 ASCII)

A	B	C	D	E	F	G	H
---	---	---	---	---	---	---	---

**Func G** Heater B (position1)

- 0) Normal
- 1) Over heating
- 2) Unusual heating time
- 3) Not available
- 4) Temperature issue

**Func H** Tip

- 0) Normal
- 1) Unusual mounting
- 2) Unusual ejecting

## 6. Response Error Code

Error Code	Description
Motor X move overtime	X-axis motor move is overtime
Motor Y move overtime	Y-axis motor move is overtime
Motor Z move overtime	Z-axis motor move is overtime
Heater heating up overtime	Heating plate heats up overtime
Heater overheated	Heating plate is overheated
Mix overtime	Mixing time exceeded
Collect overtime	Collecting time exceeded
Eject overtime	Ejecting time exceeded
Tip overtime	Tipping time exceeded
Vapor overtime	Vaporing time exceeded
TMNT FAIL	Picking up tips is failed
TEJT FAIL	Leaving tips is failed

## 6. Response Error Code

Error Code	Description
PARAM Err	Input invalid parameter value
<CMD> BUSY	Instrument is busy

## 7. Technical Support

In case of any questions, please try to contact our authorized distributor nearest to you. Taiwan Advance Nanotech Inc. provides post-sale services call number at +886-3-3167568 or via email: [service@tanbead.com](mailto:service@tanbead.com) for assistance.

Please provide this instrument serial number when you talk to our technician, that will solve the problems efficiently and answer your questions more precisely.

## 8. Cleaning and Maintenance

- Clean the device after every use. When users detect samples or reagents have been spilled, clean the device immediately to avoid damage or contamination.
- Wear gloves and appropriate personal protective equipment. If the device is used with biohazardous materials, dispose of any cleaning materials used in accordance with your institutional guidelines.
- The device may go through a run with the magnetic rods unprotected. If this happens, the magnetic rod needs to be cleaned immediately.
- To clean the magnetic rods, wipe with a soft cloth dampened with pure water. Do not use alcohol solvent.
- If the magnetic rods cannot be cleaned, please contact TANBead ([service@tanbead.com](mailto:service@tanbead.com)) for technical assistance.

## 9. Disposal

The decision whether to dispose of a potentially contaminated medical device is usually made by the owner in consultation with appropriate federal, state, and local authorities. In determining which medical devices should be discarded, the owner must assess each product's current condition and potential safety risks.

## 10. Patent

Patent List	
USA	US09616398B2
EU	EP2937136
Canada	CA2862946
Japan	JP6151735B2
Korea	KR101696517B1
China	CN104971638B
Taiwan	TWI526245B
WIPO	WO2016127292

## 11. About Manufacturer



- **Manufacturer** : Taiwan Advanced Nanotech Inc.
- **Manufacturer Address** : 6F, No. 188, Wenhe Rd. Guishan Dist., Taoyuan City 333, Taiwan
- **Manufacturer Tel** : +886-3-3167568

## 12. Advisory Notices

The Company may issue advisory notices from time to time (including safety alerts, reinforced instructions for use, or corrective recommendations) to ensure continued safe and effective use of the product.

- For the latest advisory notices, please contact Customer Service [service@tanbead.com](mailto:service@tanbead.com)
- The Company will retain records of distribution and acknowledgment of issued advisory notices and will notify distributors / end users as appropriate.

## 13. TANBead Instrument Warranty Policy

This warranty sheet covers the Nucleic Acid Extractors manufactured by Taiwan Advanced Nanotech Inc. (Hereinafter referred to as TANBead).

TANBead warrants that under normal use conditions, this product will not have any material or manufacturing defects for one year from the date of purchase. If any defect is found within the warranty period, TANBead will repair or replace the product free of charge.

This warranty policy does not apply to the following situations:

- The product is damaged due to improper use, operation, storage, maintenance, abuse, or transportation of the products.
- The product is damaged due to accident, disaster, natural phenomenon, or other force majeure factors.
- The product is modified, disassembled, reassembled, or repaired by unauthorized personnel.
- The product exceeds its expiration date or is not used according to the instructions.
- Damage or loss caused by factors beyond TANBead's control, such as sample quality, experimental conditions, or user error.

To claim the warranty, please contact our customer service department with the following information:

- Proof of purchase (invoice, receipt, etc.)
- Product name (model) and serial number
- Description of the problem and evidence of defect or malfunction (photos, test results, etc.)

We will provide you with instructions on how to return the products. Please do not return the products without our authorization. TANBead will bear the shipping cost of returning the product, but not any other expenses or losses.

This warranty policy is your only remedy and replaces any other express or implied warranties or conditions. In no event shall TANBead be liable for any indirect, special, incidental, or consequential damages, including but not limited to loss of profits, business interruption, data loss or other commercial losses.

We appreciate your business and hope that you are satisfied with our products. If you have any questions or concerns, please feel free to contact us at any time. Thank you for choosing TANBead! If you have any questions about our warranty services, please email to [service@tanbead.com](mailto:service@tanbead.com)

To register your TANBead instrument for more technical support and services, please sign up an official membership of TANBead.





Taiwan  
Advanced  
Nanotech



台灣圓點奈米技術股份有限公司  
台灣圓點核酸自動萃取儀(未滅菌)

039.L2601.X01

# Maelstrom 9610 LH








## 使用手冊

# 目錄

1.	簡介 .....	41
2.	儀器概述 .....	47
3.	安裝注意事項 .....	49
4.	指令列表 .....	50
5.	“GetErrorCode” 指令回傳對照表 .....	63
6.	錯誤訊息對照表 .....	67
7.	技術支援 .....	68
8.	清潔保養須知 .....	69
9.	儀器回收 .....	69
10.	專利資訊 .....	70
11.	製售商資訊 .....	70
12.	重要通告.....	70
13.	TANBead 儀器保固政策.....	71

# About Manual

The label on the instrument, the User Manual, and other packaging material may contain following symbols:

	型號		操作溫度限制
	序號		保持乾燥
	製造商		遠離陽光
	製造日期		歐洲授權代表
	CE標誌		體外診斷醫療器材
	使用說明書		濕度限制
	注意未滅菌		注意
	注意生物危害		生物危害
	注意高溫		
	注意機械夾手		
	廢棄物電子電機設備指令		

## 關於儀器

### 警告

- 請使用符合銷售國家標準指定電壓、電源線。
- 請勿使用損壞電源線或鬆動電源插座。
- 請勿放置液體於儀器上方，避免傾倒導致儀器損壞或故障。
- 斷開插座時請握緊儀器插頭，請勿拉扯電源線。
- 儀器清潔保養時，請關閉電源並斷開插座。
- 疑似雷擊時，請立即關閉電源並斷開插座。

## 關於儀器

### 注意事項

- 請勿撞擊、敲打儀器，造成儀器故障損壞。
- 未經製售商許可，請勿嘗試維修、改裝，造成儀器故障損壞。
- 儀器維修須經由台灣圓點奈米技術股份有限公司、授權製售商工程師操作進行。
- 儀器維修須經由台灣圓點奈米技術股份有限公司、授權製售商提供原廠零組件。
- 請依照台灣圓點奈米技術股份有限公司、授權製售商提供指示使用儀器。

# 1. 簡介

Maelstrom 9610 LH 可以整合到移液系統工作站，提高分子診斷實驗室的生產力，提供常規性操作流程轉變為無人化核酸萃取的最佳解決方案。

## 運作原理

使用磁棒收集和轉移磁珠，並使用旋轉套混合檢體及緩衝液。提昇混合效率，延伸的磁棒可以有效地收集磁珠。在細胞裂解，吸附核酸，洗滌和洗脫後獲得純化的核酸。

## 預期用途

Maelstrom 9610 LH 為純化及分離核酸之自動化操作平台，藉由磁珠從檢體中萃取、純化核酸。建議搭配使用台灣圓點奈米技術股份有限公司核酸萃取試劑相關套組，以獲得最佳的萃取純化效能。

# 1. 簡介

## 環境要求

避免減少儀器的使用壽命，請滿足以下條件使用本產品：

- 室內使用，避免溫差或濕度較大的空間
- 儀器須放置在可承載至少30公斤的平台桌面或平台上
- 操作環境要求：  
溫度: 10-40°C  
相對溼度: 40-80%
- 儲存和運輸環境要求：  
溫度: 5-50°C  
相對溼度: 20-85%
- 使用高度要求：  
低於海拔2000米

# 1. 簡介

## 安全說明及指南

- 此儀器用在具有潛在生物危害物質實驗時，需做好個人防護措施(手套、護目鏡、實驗服等)
- 在正常情況下，本產品運作時聲音不超過80分貝並造成危害，如果產生較高的分貝，請聯繫技術支援。
- 由於使用化學及生物危害物質，本設備可能存在危險。
- 使用者應遵守其所屬機構規範，處理和棄置本設備使用的所有感染性物質。
- 每次使用後需清潔儀器，若樣品或試劑打翻沾染至儀器上，請立即清潔，避免損壞或污染其他實驗樣品。
- 此儀器需與台灣圓點奈米提供之特定旋轉套一起使用，否則可能導致萃取效果不佳。
- 操作儀器前，請完整閱讀此操作說明手冊，未能閱讀、理解並遵循手冊之說明，可能導致設備損壞、操作人員受傷或實驗結果不佳。

# 1. 簡介

## 安全規範

- Maelstrom 9610 LH 設備已通過測試並符合 IEC 61010-1:2010+A1:2016 (Edition 3.1) and EN 6010-1:2010+A1:2019 《測量、控制和實驗室用電氣設備的安全要求 - 第 1 部分：一般要求》的標準。
- Maelstrom 9610 LH 設備已通過測試並符合 IEC 61010-2-101:2018 with IEC 61010-1:2010 + A1:2016 and EN 61010-2-101:2017 with EN 61010-1:2010 + A1:2019 標準《測量、控制和實驗室用電氣設備的安全要求 - 第 2-101 部分：體外診斷 (IVD) 的特殊要求》“醫用器材”。

## 電磁相容

- Maelstrom 9610 LH 設備已通過測試並符合 IEC 61326-1:2020 / EN IEC 61326-1:2021 & IEC 61326-2-6:2020 / EN IEC 61326-2-6:2021 標準。“測量、控制和實驗室用電氣設備 - EMC 要求 - 第 2-6 部分：特殊要求 - 體外診斷 (IVD) 醫療設備”。

# 1. 簡介

## 配件清單


儀器所檢附之配件，會依各國家而有所不同。



電源線



變壓器







磁座



傳輸線

## 塑膠耗材

Maelstrom 9610LH 使用專門設計的耗材，以達到最佳處理效果。使用其他類型的耗材可能會損壞儀器並影響保固。


項目	小磁棒(B)
96 旋轉套組合包	(Ø2.2)旋轉套組合包 
96 孔盤	96 孔盤 
Auto tube	16 auto tube  16 base B 

## 2.儀器概述




項目	規格
型號	Maelstrom 9610 LH
	小磁棒 (B)
磁棒	Ø2.2 mm
重量	約 25公斤
材積	542(寬) x 285(深) x 305(高)毫米
電器規格	100-240伏特, 50/60赫茲, 3.0安培 (Class I)
最大通量	96樣品/次
試劑體積	50~1,600 $\mu$ L
攪拌速度	500~3,000轉
加熱模組	1組獨立加熱盤
加熱盤	U 底
Temperate Range	RT to 100°C
Magnetic rod	>3,900 高斯

# 儀器概述



- ① 傳輸線接口
- ② LED指示燈:電源開啟/關閉
- ③ 電源開關
- ④ 電源接口

## 通訊示意圖：



### 3. 安裝注意事項

請注意本產品的重量約30公斤，從包裝箱中取出本產品進行安裝時，建議2人以上一起搬運並請特別注意安全。

**步驟 1:**

自包裝箱中取出儀器。

**步驟 2:**

移除緩衝材。

**步驟 3:**

將儀器搬運至可承載30公斤以上的桌面或是平台。

**步驟 4:**

接上電源線，請注意本產品僅適用電壓100-240伏特，錯誤的電源會導致運作異常或對儀器造成損壞。

**步驟 5:**

連接機器上USB通訊線到電腦上(Window 10 OS, RAM >8G, Intel core i5 or above)。

**步驟 6:**

打開電源，儀器將自動進行復位動作。

**步驟 7:**

操作流程會依照試劑不同而不同，請遵循試劑的使用說明書進行設定程式，一般操作流程包括：

- a. 放置旋轉套組合包到移動平台上
- b. 插取旋轉套，置換成反應盤
- c. 混合並吸取磁珠進行轉移到其他反應盤
- d. 依序置換不同的反應盤以完成萃取試驗

## 4. 指令列表

### Link Control Protocol

USB COM Port = COM50 (com0com-serial port emulators)

Baud Rate = 115200

Data Bit = 8

Stop Bit = 1

Parity Bit = None

By using \n – It prints a new line

### WELL

	Plate Position	Volume	Action	Placeholder	Placeholder
WELL	2	{Volume}	{Action}	0	0

WELL function : Set up the volume and action

{Volume} : Specified as microliter range from 50 to 1600

{Action Type} : {0} Forward

{1} Reverse

{2} Forward U/D

{3} Reverse U/D

For example :

Send to M9610LH :

"WELL 2 800 0 0 0\n" {Volume=800 $\mu$ L} · {Action=0}

M9610LH Response

"WELL ACK\n" To respond after command received

"WELL OK\n" To respond after command accomplished

Any optional field for which there is no meaningful data contains a zero (0) or a space as a placeholder.

## 4.指令列表

### PMOV

	Plate Position
PMOV	{Plate Position}

PMOV function : Set up the plate position

{Plate Position} : {1} Move to left position (only using for 48 series)  
{2} Move to right position  
{3} Move to exchange position  
{4} Distance from initial point 2 mm

For example :

Send to M9610LH :

"PMOV 2\n"                      {Plate Position = 2}

M9610LH Response :

"PMOV ACK\n"

To respond after command received.

"GO\_FRONT OK\n"

To respond after shield plate moving to front limit position.

"PMOV OK\n"

To respond after command accomplished.

## 4.指令列表

### PMIX

	Plate Position	Spin Time	Spin Speed
PMIX	2	{Spin Time}	{Spin Speed}

PMIX function : Set up spin time and speed

{Spin Time} : Specified as minutes from 0.1 to 60

{Spin Speed} : Specified as RPM from 500 to 3,000

For example :

Send to M9610LH :

"PMIX 2 0.5 3000\n"

{Spin Time = 0.5(min)}

{Spin Speed = 3000(rpm)}

M9610LH Response :

"PMIX ACK\n"

To respond after command received.

"GO\_REAR OK\r\n"

To respond after shield plate moving to rear limit position.

"PMIX GO\n"

To respond after moving to the position and start the real mix action.

"GO\_FRONT OK\r\n"

The mix time would start to count down.

To respond after shield plate moving to front limit position.

"GO\_REAR OK\r\n"

To respond after shield plate moving to rear limit position.

"PMIX OK\n"

To respond after command accomplished.

## 4.指令列表

### PCOL

	Plate Position	Beads Collection
PCOL	2	{Collect Time}

PCOL function : Set up the collective time

{Collect Time} : Specified as minutes from 0.1 to 60

For example :

Send to M9610LH :

"PCOL 2 0.5\n"                    {Collect time = 0.5(min)}

M9610LH Response :

"PCOL ACK\n"

To respond after command received.

"GO\_REAR OK\n"

To respond after shield plate moving to rear limit position

"PCOL POS1\n"

To respond after moving to the position 1 and start the real collect action.

"PCOL POS2\n"

The collect time would start to count down.  
To respond after moving to the position 2 and start the real collect action.

"PCOL POS3\n"

The collect time would start to count down.  
To respond after moving to the position 3 and start the real collect action.

"PCOL POS4\n"

The collect time would start to count down.  
To respond after moving to the position 4 and start the real collect action.

"PCOL POS5\n"

The collect time would start to count down.  
To respond after moving to the position 5 and start the real collect action.

"GO\_FRONT OK\n"

The collect time would start to count down.  
To respond after shield plate moving to front limit position.

"GO\_REAR OK\n"

To respond after shield plate moving to rear limit position.

"PCOL OK\n"

To respond after command accomplished.

## 4.指令列表

### PVAPOR

	Plate Position	Vapor time
PVAPOR	2	{Vapor Time}

PVAPOR Function : Set up the vapor time

{Vapor Time} : Specified as minutes from 0.1 to 60

For example :

Send to M9610LH :

"PVAPOR 2 0.5\n"

{Vapor time=0.5(min)}

M9610LH Response:

"PVAPOR ACK\n"

To respond after command received.

"GO\_FRONT OK\r\n"

To respond after shield plate moving to front limit position and start the real vapor action. The vapor time would start to count down.

"GO\_REAR OK\r\n"

To respond after shield plate moving to rear limit position.

"GO\_FRONT OK\r\n"

To respond after shield plate moving to front limit position.

"PVAPOR OK\n"

To respond after command accomplished.

## 4.指令列表

### TMNT

	Plate Position
TMNT	2

TMNT function : Pick up tips  
For example:

Send to M9610LH :  
"TMNT 2\n"

M9610LH Response :

"TMNT ACK\n"

"GO\_REAR OK\r\n"

"GO\_FRONT OK\r\n"

"TMNT OK\n"

To respond after command received.

To respond after shield plate moving to rear limit position.

To respond after shield plate moving to front limit position.

To respond after command accomplished.

## 4.指令列表

### TEJT

	Plate Position
TEJT	2

TEJT function : Leave tips

For example :

Send to M9610LH :

"TEJT 2\n"

M9610LH Response :

"TEJT ACK\n"

"GO\_REAR OK\r\n"

"GO\_FRONT OK\r\n"

"TEJT OK\n"

To respond after command received.

To respond after shield plate moving to rear limit position.

To respond after shield plate moving to front limit position.

To respond after command accomplished.

## 4.指令列表

### TMPE

	Placeholder	Temperature
TMPE	0	{Temperature}

TMPE function : Set up the temperature

{Temperature} : Specified as degrees Celsius from RT to 130

For example :

Send to M9610LH :

"TMPE 0 40\n" {Temperature=40 (deg C)}

M9610LH Response :

"TMPE ACK\n" To respond after receiving command and start the heating process if necessary.

"TMPE OK\n" To respond after temperature reaches the target or after turning off the heating process.

Caution: Please turn off TMPE after heating process.

Send to M9610LH :

"TMPE 0 0\n" {Temperature = 0 (turn Off)}

M9610LH Response :

"TMPE ACK\n" To respond after receiving command and start the heating process if necessary.

"TMPE OK\n" To respond after temperature reaches the target or after turning off the heating process.

## 4.指令列表

### PAUS

	Pause or Resume
PAUS	{Pause or Resume}

PAUSE function : PAUSE is only using during spinning period

{Pause or Resume} : {0} Resume  
{1} Pause

For example :

Send to M9610LH :  
"PMIX 2 0.5 3000\n"

M9610LH Response :

"PMIX ACK\n"

"GO\_REAR OK\r\n"

"PMIX GO\n"

To respond after command received.

To respond after shield plate moving to rear limit position

To respond after moving to the position and start the real mix action. The mix time would start to count down.

Send to M9610LH :

"PAUS 1\n"

pause mix command

M9610LH Response :

"PAUS ACK\n"

"PAUS OK\n"

"GO\_FRONT OK\r\n"

To respond after command received. The mix time would pause to count down.

To respond after command accomplished.

To respond after shield plate moving to front limit position

## 4.指令列表

### PAUS

Send to M9610LH :

"PAUS 0\n" resume mix command

M9610LH Response :

"PAUS ACK\n" To respond after command received.

"PAUS OK\n" To respond after command accomplished.

"GO\_REAR OK\n" To respond after shield plate moving to rear limit position.

"PMIX GO\n" To respond after moving to the position and start the real mix action. The mix time would resume to count down.

"PMIX OK\n" To respond after command accomplished.

### STOP

STOP function: The instrument will be reset and disconnected. And this function is only used for spinning period and idle.

For example :

Send to M9610LH :

"STOP\n"

M9610LH Response :

"STOP ACK\n" To respond after command received.

"STOP OK\n" To respond after command accomplished.

p.s. When the instrument is executing a STOP order, the stop process needs to be accomplished. Otherwise, the operator can't give a command.

### ORIGIN

ORIGIN function: Three axis (X, Y and Z) moving to basis point

For example :

Send to M9610LH :

"ORIGIN\n"

M9610LH Response :

"ORIGIN ACK\n" To respond after command received.

"ORIGIN OK\n" To respond after command accomplished.

## 4.指令列表

### Read Temp

ReadTemp function : Reading temperature

For example :

Send to M9610LH :

"ReadTemp\n"

M9610LH Response :

"ReadTemp 0.00 19.79\n"

To respond after command received.

"ReadTemp {Placeholder (48LH only)} {Temperature}"

### GetFWVer

GetFWVer function: Getting firmware version

For example :

Send to M9610LH :

"GetFWVer\n"

M9610LH Response :

"GetFWVer V1.0.0.2T1\_9610LH\n"

To respond after command received.

"GetFWVer {FW Version}"

## 4.指令列表

### TMPR

	Switch Temperature Sensor
TMPR	{Temperature Sensor Type}

TMPR function: Control the temperature sensor

{Switch Temperature Sensor} :    {0} Close  
  {1} Open

For example :

Send to M9610LH :  
"TMPR 1\n"                    Open TMPR

M9610LH Response :  
"TMPR ACK\n"                To respond after command received.

Send to M9610LH :  
"ReadTemp\n"                Read temperature command

M9610LH Response :  
"ReadTemp 478.54 48.18\n"  
                                  To respond after command received.  
                                  {Temperature = 48.18}

Send to M9610LH :  
"ReadTemp\n"                Read temperature command

M9610LH Response :  
"ReadTemp 478.55 48.18\n"  
                                  To respond after command received.  
                                  {Temperature = 48.18}

Send to M9610LH :  
"ReadTemp\n"                Read temperature command

M9610LH Response :  
"ReadTemp 478.55 44.47\n"  
                                  To respond after command received.  
                                  {Temperature = 44.47}

## 4.指令列表

### TMPR

Send to M9610LH :  
"ReadTemp\n"                      Read temperature command

M9610LH Response :  
"ReadTemp 478.55 44.40\n"  
To respond after command received.  
{Temperature= 44.40}

Send to M9610LH :  
"TMPR 0\n"                      Close TMPR

M9610LH Response :  
"TMPR ACK\n"                      To respond after command received.  
"TMPR OK\n"                      To respond after command accomplished.

### GetErrorCode

GetErrorCode function: Getting error code

For example :  
Send to M9610LH :  
"GetErrorCode\n"

M9610LH Response :  
"GetErrorCode 1111111 00000030\n"  
To respond after command received.  
"GetErrorCode {state code(7 ASCII)} {error code(8 ASCII)}"

## 5. “GetErrorCode” 指令回傳對照表

1<sup>st</sup> argument: State code (7 ASCII)

1	2	3	4	5	6	7
---	---	---	---	---	---	---

### Func 1 Plate state

- 1) Basis point
- 2) Remove position
- 3) Left position
- 4) Right position
- 5) Moving
- 6) Abnormal

### Func 2 Mixing state

- 1) Normal
- 2) Spinning
- 3) Pause
- 4) Abnormal

### Func 3 Collection state

- 1) Normal
- 2) Elution
- 3) Collection
- 4) Abnormal

### Func 4 Heater state A (position2)

- 1) Stop heating
- 2) Reach setting temperature
- 3) Heating up
- 4) Placeholder
- 5) Abnormal

### Func 5 Heater state B (position1) Func 6 Tip state

- 1) Stop heating
- 2) Reach setting temperature
- 3) Heating up
- 4) Placeholder
- 5) Abnormal

- 1) Normal
- 2) Mounting
- 3) Ejecting
- 4) Abnormal

## 5. "GetErrorCode" 指令回傳對照表

1<sup>st</sup> argument: State code (7 ASCII)

1	2	3	4	5	6	7
---	---	---	---	---	---	---

**Func 7** Vapor state

- 1) Normal
- 2) Vapor
- 3) Abnormal

## 5. "GetErrorCode" 指令回傳對照表

2<sup>nd</sup> argument: Error code (8 ASCII)

A	B	C	D	E	F	G	H
---	---	---	---	---	---	---	---

**Func A** X-axis motor

- 0) Normal
- 1) Placeholder
- 2) Limit deviate
- 3) Calibration failed

**Func B** Y-axis motor

- 0) Normal
- 1) Placeholder
- 2) Placeholder
- 3) Calibration failed

**Func C** Z-axis motor

- 0) Normal
- 1) Placeholder
- 2) Placeholder
- 3) Calibration failed

**Func D** Slew drive motor

- 0) Normal
- 1) Unusual rotational speed

**Func E** Drip-proof motor

- 0) Normal
- 1) Unusual front position
- 2) Unusual rear position

**Func F** Heater A (position2)

- 0) Normal
- 1) Over heating
- 2) Unusual heating time
- 3) Temperature issue

## 5. “GetErrorCode” 指令回傳對照表

2<sup>nd</sup> argument: Error code (8 ASCII)

A	B	C	D	E	F	G	H
---	---	---	---	---	---	---	---

**Func G** Heater B (position1)

- 0) Normal
- 1) Over heating
- 2) Unusual heating time
- 3) Not available
- 4) Temperature issue

**Func H** Tip

- 0) Normal
- 1) Unusual mounting
- 2) Unusual ejecting

## 6. 錯誤訊息對照表

錯誤代碼	解析
Motor X move overtime	X軸馬達移動逾時
Motor Y move overtime	Y軸馬達移動逾時
Motor Z move overtime	Z軸馬達移動逾時
Heater heating up overtime	加熱盤加熱逾時
Heater overheated	加熱盤過熱
Mix overtime	攪拌逾時
Collect overtime	磁吸逾時
Eject overtime	退套逾時
Tip overtime	插套逾時
Vapor overtime	風乾逾時
TMNT FAIL	插套異常
TEJT FAIL	退套異常

## 6. 錯誤訊息對照表

錯誤代碼	解析
PARAM Err	輸入無效參數值
<CMD> BUSY	設備忙碌

## 7. 技術支援

本產品提供完整的售後服務及技術支援。若使用上有任何問題，請嘗試聯繫離您最近的授權經銷商或原廠製售商。

電話: +886-3-3167568

信箱: [service@tanbead.com](mailto:service@tanbead.com)

請在與我們的技術人員交談時提供此儀器序列號，將更有效地解決、回答問題。

## 8. 清潔保養須知

- 每次實驗後，當用戶檢測到樣品或試劑溢出時和請立即清潔設備，以免損壞或污染。
- 請戴上手套和適當的個人防護設備清潔儀器，並按照您的機構指南處理使用後的清潔用品。
- 儀器可能會在磁棒不受保護的情況下運行，發生這種情況，請立即清潔磁棒。
- 請使用純淨水和軟布擦拭清潔磁棒，請勿使用酒精溶劑。
- 磁棒無法清潔時，請聯繫TANBead（電話：+886-3-3167568）尋求技術支援。

## 9. 儀器回收

適當處置可能受污染的醫療器械的決定通常由所有者與適當的聯邦、州和地方當局協商後做出。在確定應丟棄哪些醫療設備時，所有者必須評估每個產品的當前狀況和潛在的安全風險。

## 10. 專利資訊

### 專利列表

美國	US09616398B2
歐洲	EP2937136
加拿大	CA2862946
日本	JP6151735B2
韓國	KR101696517B1
中國	CN104971638B
台灣	TWI526245B
世界智慧財產權組織	WO2016127292

## 11. 製售商資訊



- 製售商名稱：台灣圓點奈米技術股份有限公司
- 製售商地址：333桃園市龜山區文禾路188號6樓
- 製售商電話：+886-3-3167568

## 12. 重要通告資訊

本公司可能會不定期發佈重要通告（包含安全提示、操作強化說明或產品修正建議），以保障使用安全與產品效能。

- 最新通告請聯絡本公司客服窗口：[service@tanbead.com](mailto:service@tanbead.com)
- 本公司對已發佈之通告將保留分發與簽收紀錄，並在必要時通知代理商與使用單位採取相應行動。

## 13. TANBead 儀器保固政策

本保固條款適用於台灣圓點奈米技術股份有限公司（以下簡稱 TANBead）所製造之核酸萃取儀。

TANBead 保證於正常使用條件下，本產品自購買日起一年內不會產生材料或製造上的缺陷。如於保固期間內發現產品有任何瑕疵，TANBead 將提供免費維修或更換服務。

以下情況不適用於本保固政策：

- 因使用、操作、儲存、保養不當、誤用或運輸過程中造成之損壞
- 因意外、災害、天然現象或其他不可抗力因素造成之損壞
- 經未授權人員改裝、拆解、重組或維修者
- 產品超過其使用期限或未依使用說明操作
- 因非 TANBead 可控因素（如樣本品質、實驗條件或操作錯誤）所導致之損壞或損失

提出保固申請時，請提供以下資訊：

- 購買憑證（發票、收據等）
- 產品名稱（型號）與序號
- 問題描述及相關證明（如照片、測試結果等）

我們將提供產品退回流程指引。請勿在未獲授權情況下直接退回產品。TANBead 將負擔產品退回之運費，但不承擔其他相關費用或損失。

本保固政策為您唯一之補償條款，並取代其他明示或暗示之保證條件。TANBead 不對任何間接、特殊、附帶或衍生性損害負責，包括但不限於利潤損失、業務中斷、資料遺失或其他商業損失。

感謝您選擇 TANBead，我們致力於提供讓您滿意的產品與服務。若有任何疑問，歡迎隨時聯絡我們（[service@tanbead.com](mailto:service@tanbead.com)）。

如欲獲得更多技術支援與服務，請註冊加入 TANBead 官方會員。





Taiwan  
Advanced  
Nanotech

**Taiwan Advanced Nanotech Inc.** 

---

6F., No. 188, Wenhe Rd., Guishan Dist., Taoyuan City 333, Taiwan (R.O.C.)

Tel : +886-3-316-7568 Fax: +886-3-317-3369

[www.tanbead.com](http://www.tanbead.com)