



# TANBead® Nucleic Acid Extraction Kit

## Viral Auto Tube

(For use with the SLA-16 / 32 / E13200 series)



635S46

(For Professional Use Only) V3

### 1. Intended Purpose

The TANBead® Nucleic Acid Extraction Kit is a nucleic acid purification kit based on magnetic bead technology by using with corresponding. TANBead® Nucleic Acid Extraction Kit is suitable for isolating a variety of viral nucleic acids such as Enterovirus, Japanese encephalitis virus, Dengue virus, Avian Influenza virus, and Epstein-Barr virus. Serum specimens are processed through a series of automatic extraction steps and finally the high-quality nucleic acids can be applied directly to the following qualitative and quantitative assays. The kit is intended for use by technicians, physicians, and biologists with well-trained in molecular biological techniques, the techniques of magnetic bead purification and *in vitro* diagnostic procedures. Any diagnostic results generated by using the sample preparation procedure in conjunction with any downstream diagnostic assay should be interpreted related to other clinical or laboratory findings. The kit is not limited to any specific disorder, condition, or other additional accompanying diagnostics. It is applicable for all population.

### 2. The basic principle

The silicon dioxide layer coated on the magnetic beads can adsorb the negatively charged molecules to purify nucleic acids from samples.

### 3. Specification

|                    |                                |
|--------------------|--------------------------------|
| Starting Materials | 200 µL serum or PBS suspension |
| Elution Volume     | 70~100 µL                      |

### 4. Component Supplied with the Kit

|                |         |                                |
|----------------|---------|--------------------------------|
| Auto Tube      | 6 trays | Auto Tube with reagent buffers |
| Proteinase K   | 1.0 mL  | Proteinase K                   |
| Elution Buffer | 20 mL   | Nuclease-Free Water            |
| Strip          | 12      | 8-channel strip                |
| Protocol       | 1       | Instruction guide for user     |

### 5. Auto Tube Content

| Well   | Buffer           | Volume (µL) |
|--------|------------------|-------------|
| 1 / 7  | Lysis Buffer     | 700         |
| 2 / 8  | Washing Buffer 1 | 800         |
| 3 / 9  | Magnetic Beads   | 800         |
| 4 / 10 | Washing Buffer 2 | 800         |
| 5 / 11 | Washing Buffer 2 | 800         |
| 6 / 12 | Elution Buffer   | 100         |

### 6. Kit Storage and Shelf Life

- Components under room temperature (15~35°C) can be stored until the expiration date labeled on the box.
- The Proteinase K was transported at room temperature. When received, please store at 4°C.

### 7. Precautions

- It can only be used for *in vitro* diagnostic.
- Avoid using expired reagents.
- When the temperature is below 20°C, place the Auto Plates / Auto Tubes in an oven (preheated 42~60°C) 5 to 10 minutes.
- Avoid vigorous shaking, in order to avoid excessive formation of foam.
- Carefully remove aluminum foil to avoid splashing.
- Do not expose the opened reagents or Auto Plates / Auto Tubes to air. The evaporation would lead to pH change, or effect on the extraction effectiveness.
- Please check the integrity of the Auto Plates / Auto Tubes and remember to insert the strips into the appropriate position of the suitable instrument before operating them.
- Please wear a mask and disposable gloves when handling.
- Use sterile consumables to avoid nuclease contamination.
- Reagent solution contains guanidine salt, avoid using bleach containing detergent.

- Avoid eyes, skin, and clothing contact with reagents. In case of any contact, flush with flowing water.
- If any serious incident occurs, please report to the manufacturer and the competent authority of the member state in which the user and / or the patient is established.

### 8. Materials required, Not Supplied

- TANBead® Nucleic Acid Extraction System  
Model: SLA-16 / 32 / E13200 series (non-sterile)
- Disposable gloves
- Scissors, utility knives
- Micropipette, disposable tips (10 µL / 200 µL / 1000 µL)
- 1.5 mL microcentrifuge tube
- 15 mL / 50 mL conical tube

### 9. Sample Collection, Transportation, and Storage

#### ■ Sample collection and storage

- Serum, whole blood
  - Serum specimens must be obtained from serum collection tubes, whole blood specimens must be obtained from sodium citrate or EDTA collection tubes.
- Specimen storage
  - Fresh whole blood specimens can be stored at room temperature for 6 hours.
  - After centrifugation, the serum sample can be stored at
    - Room temperature for 24 hours.
    - 2~8°C up to 7 days.
    - 20°C for long-term preservation.

#### ■ Specimen transportation

Transportation of whole blood, serum specimens should be followed by specific pathogen transportation-related laws. The whole blood sample should be kept between 2~25°C during transportation and within 6 hours for separated serum. Serum samples can be transported between 2~8°C or by freezing.

### 10. Nucleic Acids Extraction Protocol

Before operating, turn on the warm-up system of TANBead® Nucleic Acid Extractor, if it is equipped with temperature controller, please setting at **45°C**.

- Carefully remove the aluminum foil on the Auto Tubes.
- Pipet **200 µL serum or PBS suspension** and **10 µL Proteinase K** into column **#1 / #7** of Auto Tube.

**Note: The volume ratio of mixture and lysis buffer is about 200 µL : 700 µL. If it is changed, it might be affected the performance.**

- Push Assembled Auto Tubes completely to the bottom of the plate rack. Make sure that the chamfer of the plate is at the lower left.
- Push strips completely to the bottom of strip rack frame.
- Close the door panel.
- Select the program "**VIRUS-W4-AUTO**". The parameters are given in following section.
- Carefully remove the Auto Tubes when the program is finished.
- Use micropipette to transfer the purified nucleic acid from column **#6 / #12** to a clean tube.
- Discard used Auto Tubes and strips into the waste recycling bin.

## 11. Program

### ■ SLA-16 / 32 series

| Program Name: VIRUS-W4-AUTO |      |            |             |     | Model: SLA-16 / 32 series |             |       |           |
|-----------------------------|------|------------|-------------|-----|---------------------------|-------------|-------|-----------|
| Step                        | Well | Mixing (M) | Collect (S) | Rod | Mixing speed              | Volume (μL) | Pause | Vapor (M) |
| 1                           | 3    | 1          | 60          | On  | Medium                    | 800         | Off   | 0         |
| 2                           | 2    | 1          | 0           | Off | Medium                    | 800         | Off   | 0         |
| 3                           | 1    | 20         | 0           | Off | Low                       | 900         | Off   | 0         |
| 4                           | 2    | 0          | 60          | On  | Medium                    | 800         | Off   | 0         |
| 5                           | 1    | 10         | 60          | On  | Medium                    | 900         | Off   | 0         |
| 6                           | 2    | 2          | 60          | On  | Medium                    | 800         | Off   | 0         |
| 7                           | 3    | 2          | 60          | On  | Medium                    | 800         | Off   | 0         |
| 8                           | 4    | 2          | 60          | On  | Medium                    | 800         | Off   | 0         |
| 9                           | 5    | 2          | 60          | On  | Medium                    | 800         | Off   | 10        |
| 10                          | 6    | 5          | 120         | On  | Medium                    | 150         | Off   | 0         |
| 11                          | 5    | 1          | 0           | Off | Medium                    | 800         | Off   | 0         |
| 12                          | 0    | 0          | 0           | Off | Medium                    | 0           | Off   | 0         |

### ■ SLA-E13200 series

| Program Name: VIRUS-W4-AUTO |      |           |            |             |     | Model: SLA-E13200 series |             |       |           |
|-----------------------------|------|-----------|------------|-------------|-----|--------------------------|-------------|-------|-----------|
| Step                        | Well | Temp (°C) | Mixing (M) | Collect (S) | Rod | Mixing speed             | Volume (μL) | Pause | Vapor (M) |
| 1                           | 3    | 45        | 1          | 60          | On  | Medium                   | 800         | Off   | 0         |
| 2                           | 2    | 45        | 1          | 0           | Off | Medium                   | 800         | Off   | 0         |
| 3                           | 1    | 45        | 20         | 0           | Off | Low                      | 900         | Off   | 0         |
| 4                           | 2    | 45        | 0          | 60          | On  | Medium                   | 800         | Off   | 0         |
| 5                           | 1    | 45        | 10         | 60          | On  | Medium                   | 900         | Off   | 0         |
| 6                           | 2    | 45        | 2          | 60          | On  | Medium                   | 800         | Off   | 0         |
| 7                           | 3    | 45        | 2          | 60          | On  | Medium                   | 800         | Off   | 0         |
| 8                           | 4    | 45        | 2          | 60          | On  | Medium                   | 800         | Off   | 0         |
| 9                           | 5    | 45        | 2          | 60          | On  | Medium                   | 800         | Off   | 10        |
| 10                          | 6    | 45        | 5          | 120         | On  | Medium                   | 150         | Off   | 0         |
| 11                          | 5    | N / A     | 1          | 0           | Off | Medium                   | 800         | Off   | 0         |
| 12                          | 0    | N / A     | 0          | 0           | Off | Medium                   | 0           | Off   | 0         |

## 12. Result

Nucleic acid product purified by TANBead® nucleic acid extraction kit can perform qualitative / quantitative analysis of specific genes by PCR, RT-PCR, Q-PCR or Qrt-PCR. Please refer to the molecular diagnostic kit manual.

## 13. Reagent performance

### ■ Repeatability

A five-day reproducibility test was carried out with 100 IU / mL of EBV serum samples for 5 consecutive days with the same reagent kit by different operators. The coefficient of variation of nucleic acid extraction concentration is less than 5%.

### ■ Reproducibility

A five-day reproducibility test was carried out with 100 IU / mL of EBV serum samples for 5 consecutive days with the same reagent kit by different operators. The coefficient of variation of nucleic acid extraction concentration is less than 5%.

### ■ Detection limit of HCV virus: ≥100 IU / mL

### ■ Interfering substance

According to preclinical tests, the performance of extraction kit will not be affected by EDTA, Li-Heparin, Sodium Citrate, D-Glucose, Hemoglobin, lipoprotein and triglyceride in samples.

## ■ The stability of extracted RNA

| Storage Conditions | RNA stability |
|--------------------|---------------|
| -80°C              | Over 90 days  |
| -20°C              | 28 days       |
| 4°C                | 14 days       |
| 25°C               | 2 days        |
| Freeze-thaw        | 5 times       |

## 14. Explanation of Symbols

|  |                     |  |   |
|--|---------------------|--|---|
|  | Manufacturer        |  | Consult instructions for use              |
|  | Temperature limit   |  | Contains sufficient for test              |
|  | CE mark             |  | In vitro diagnostic medical use           |
|  | Catalogue number    |  | Caution                                   |
|  | Batch code          |  | Non-sterile                               |
|  | Do not re-use       |  | Protect from heat and radioactive sources |
|  | Date of manufacture |  | Use-by date                               |

EC REP

mdi Europa GmbH, Langenhagener Str. 71, 30855 Langenhagen, Germany

## 15. Post-market surveillance conclusion

After a risk assessment and clinical evaluation assessment, when weighing the benefits of medical device, patients, and the risks associated with the use of the device, the risk is acceptable. The post-market surveillance report shows that no death or serious adverse events occurred.