



# TANBead® Nucleic Acid Extraction Kit

## Rice DNA Auto Tube

(For use with the SLA-16 / 32 / E13200 series)

RUO

619S45

(For Research Use Only) V3

### 1. Intended Use

This product is designed for isolating nucleic acid from various samples, which can be performed by using TANBead® Nucleic Acid Extractor and is intended for research use only.

### 2. Purpose

TANBead® Nucleic Acid Extraction Kit (619S45) with TANBead® Nucleic Acid Extractor (SLA-16/32, SLA-E132 Series) simplifies the process of isolating rice DNA from grinding grain and cell lysis to DNA purification completely by automated system. In addition, because TANBead has its own nano-beads for nucleic acid binding with high efficiency, it takes only one single grain of rice for nucleic acid extraction and is ideal for applications such as germplasm identification and laboratories of academic research. TANBead® Nucleic Acid Extractor can handle 32 rice samples simultaneously and the automated process can save labor and time by enhancing the extraction efficiency, consistency and reproducibility.

### 3. The basic principle

The silicon dioxide layer coated on the magnetic beads can adsorb the negatively charged molecules to purify nucleic acids from samples.

### 4. Specification

|                    |                          |
|--------------------|--------------------------|
| Starting Materials | One single grain of rice |
| Elution Volume     | 90~130 µL                |

### 5. Component Supplied with the Kit



|                |            |  |
|----------------|------------|--|
| Auto Tube      | 4 trays    | Auto Tube with reagent buffers           |
| Lysis Buffer 1 | 25 mL x 1  | Tris buffer, surfactants, pH 8.0         |
| Lysis Buffer 2 | 80 mL x 1  | Guanidine salt, Tris buffer, surfactants |
| Elution Buffer | 20 mL x 1  | Nuclease-Free Water                      |
| Proteinase K   | 1.0 ml x 1 | Proteinase K                             |
| Base           | 2          | A rack for 8 Auto Tubes                  |
| Strip          | 12         | 8-channel strip                          |
| Protocol       | 1          | Instruction guide for user               |

### 6. Auto Tube Content

| Well   | Buffer           | Volume (µL) |
|--------|------------------|-------------|
| 1 / 7  | Elution Buffer   | 130         |
| 2 / 8  | Washing Buffer 3 | 800         |
| 3 / 9  | Magnetic Beads   | 800         |
| 4 / 10 | Washing Buffer 1 | 800         |
| 5 / 11 | -                | -           |
| 6 / 12 | -                | -           |

### 7. Kit Storage and Shelf Life

- Components under room temperature (15~35°C) can be stored until the expiration date labeled on the box.
- The proteinase K is transported at room temperature. Upon received, please store proteinase K at 2~8°C.

### 8. Precautions

- For research use only.
- Avoid using expired reagents.
- When the temperature is below 20°C, place the Auto Plates / Auto Tubes in an oven (preheated 42~60°C) 5 to 10 minutes.
- Avoid vigorous shaking, in order to avoid excessive formation of foam.
- Carefully remove aluminum foil to avoid splashing.
- Do not expose the opened reagents or Auto Plates / Auto Tubes to air. The evaporation would lead to pH change, or effect on the extraction effectiveness.
- Please check the integrity of the Auto Plates / Auto Tubes and remember to insert the strips into the appropriate position of the suitable instrument before operating them.
- Please wear a mask and disposable gloves when handling.
- Use sterile consumables to avoid nuclease contamination.

- Reagent solution contains guanidine salt, avoid using bleach containing detergent.
- Avoid eyes, skin, and clothing contact with reagents. In case of any contact, flush with flowing water.
- If any serious incident occurs, please report to the manufacturer and the competent authority of the member state in which the user and / or the patient is established.

### 9. Materials required, Not Supplied

- TANBead® Nucleic Acid Extraction System  
Model: SLA-16 / 32 / E13200 series (non-sterile)
- Disposable gloves
- Scissors, utility knives
- Micropipette, disposable tips (10 µL / 200 µL / 1000 µL)
- 1.5 mL microcentrifuge tube
- 15 mL / 50 mL conical tube

### 10. Nucleic Acids Extraction Protocol

Before operating, turn on the warm-up system of TANBead® Nucleic Acid Extractor, if it is equipped with temperature controller, please setting at **70°C**.

- Preparing the Assembled Auto Tubes by inserting Auto Tubes into the Base completely.
- Carefully remove the aluminum foil on the Auto Tubes.
- Use micropipette to load **400 µL** Lysis Buffer 1 and **20 µL** Proteinase K into column **#6 / #12** of Auto Tube.
- Put one naked rice seed into column **#6 / #12** of Auto Tube.
- Push Assembled Auto Tubes combined with conducting plate which is attached to column **#6 / #12** completely to the bottom of plate rack completely to the bottom of the plate rack. Make sure that the chamfer of the plate is at the lower left.
- Push strips completely to the bottom of strip rack frame.
- Close the door panel.
- Select the program "**DNA RICE AUTO**". The parameters are given in following section.
- The program will stop temporarily and the buzzer alarm after one hour.
- If Nucleic Acid Extractor is equipped with temp. controller, please set at **45°C**.
- Carefully pull out Assembled Auto Tubes and pay attention to high temperature of conducting plate.
- Use micropipette to transfer half of the lysate from column **#6 / #12** to **#5 / #11** and load **750 µL** Lysis Buffer 2 into **#5 / #6** and **#11 / #12**.
- Push Assembled Auto Tubes completely to the bottom of plate rack, make sure that chamfer of the plate is at the lower left and press "Pause" to go on program.
- Carefully remove the Auto Tubes when the program is finished.
- Use micropipette to transfer the purified nucleic acids from well **#6 / #12** to a clean tube.
- Discard used Auto Tubes and strips into the waste recycling bin.

## 11. Program

### ■ SLA-16 / 32 series

| Program Name: DNA-RICE-AUTO |      |            |             |     | Model: SLA-16 / 32 series |             |       |           |
|-----------------------------|------|------------|-------------|-----|---------------------------|-------------|-------|-----------|
| Step                        | Well | Mixing (M) | Collect (S) | Rod | Mixing speed              | Volume (μL) | Pause | Vapor (M) |
| 1                           | 6    | 20         | 0           | Off | Medium                    | 600         | Off   | 0         |
| 2~26                        | 6    | 1          | 0           | Off | Medium                    | 600         | Off   | 0         |
| 27                          | 6    | 10         | 0           | Off | Medium                    | 600         | Off   | 0         |
| 28                          | 6    | 5          | 0           | Off | Medium                    | 600         | On    | 0         |
| 29                          | 6    | 5          | 0           | Off | Medium                    | 900         | Off   | 0         |
| 30                          | 5    | 5          | 0           | Off | Medium                    | 900         | Off   | 0         |
| 31                          | 3    | 1          | 60          | On  | Medium                    | 800         | Off   | 0         |
| 32                          | 4    | 1          | 60          | On  | Medium                    | 800         | Off   | 0         |
| 33                          | 6    | 5          | 100         | On  | Medium                    | 900         | Off   | 0         |
| 34                          | 5    | 5          | 100         | On  | Medium                    | 900         | Off   | 0         |
| 35                          | 4    | 3          | 100         | On  | Medium                    | 800         | Off   | 0         |
| 36                          | 3    | 2          | 60          | On  | Medium                    | 800         | Off   | 0         |
| 37                          | 2    | 2          | 60          | On  | Medium                    | 800         | Off   | 10        |
| 38                          | 1    | 5          | 100         | On  | Medium                    | 150         | Off   | 0         |
| 39                          | 2    | 1          | 0           | Off | Medium                    | 800         | Off   | 0         |
| 40                          | 0    | 0          | 0           | Off | Medium                    | 0           | Off   | 0         |

### ■ SLA-E13200 series














| Program Name: DNA-RICE-AUTO |      |           |            |             |     | Model: SLA-E13200 series |             |       |           |
|-----------------------------|------|-----------|------------|-------------|-----|--------------------------|-------------|-------|-----------|
| Step                        | Well | Temp (°C) | Mixing (M) | Collect (S) | Rod | Mixing speed             | Volume (μL) | Pause | Vapor (M) |
| 1                           | 6    | 70        | 20         | 0           | Off | Medium                   | 600         | Off   | 0         |
| 2~26                        | 6    | 70        | 1          | 0           | Off | Medium                   | 600         | Off   | 0         |
| 27                          | 6    | 70        | 10         | 0           | Off | Medium                   | 600         | Off   | 0         |
| 28                          | 6    | 45        | 5          | 0           | Off | Medium                   | 600         | On    | 0         |
| 29                          | 6    | 45        | 5          | 0           | Off | Medium                   | 900         | Off   | 0         |
| 30                          | 5    | 45        | 5          | 0           | Off | Medium                   | 900         | Off   | 0         |
| 31                          | 3    | 45        | 1          | 60          | On  | Medium                   | 800         | Off   | 0         |
| 32                          | 4    | 45        | 1          | 60          | On  | Medium                   | 800         | Off   | 0         |
| 33                          | 6    | 45        | 5          | 100         | On  | Medium                   | 900         | Off   | 0         |
| 34                          | 5    | 45        | 5          | 100         | On  | Medium                   | 900         | Off   | 0         |
| 35                          | 4    | 45        | 3          | 100         | On  | Medium                   | 800         | Off   | 0         |
| 36                          | 3    | 45        | 2          | 60          | On  | Medium                   | 800         | Off   | 0         |
| 37                          | 2    | 45        | 2          | 60          | On  | Medium                   | 800         | Off   | 10        |
| 38                          | 1    | 45        | 5          | 100         | On  | Medium                   | 150         | Off   | 0         |
| 39                          | 2    | N / A     | 1          | 0           | Off | Medium                   | 800         | Off   | 0         |
| 40                          | 0    | N / A     | 0          | 0           | Off | Medium                   | 0           | Off   | 0         |

## 12. Reagent performance

### ■ The stability of extracted DNA

| Storage Conditions | DNA stability |
|--------------------|---------------|
| -80°C              | Over 90 days  |
| -20°C              | 28 days       |
| 4°C                | 14 days       |
| 25°C               | 2 days        |
| Freeze-thaw        | 10 times      |

## 13. Explanation of Symbols

|   |                       |   |                              |
|---|-----------------------|---|------------------------------|
|  | Manufacturer          |  | Consult instructions for use |
|  | Temperature limit     |  | Contains sufficient for test |
|  | Catalogue number      |  | Caution                      |
|  | Batch code            |  | Non-sterile                  |
|  | Do not re-use         |  | Keep away from sunlight      |
|  | Date of manufacture   |  | Use-by date                  |
|  | For research use only |   |                              |

