



TANBEAD

TANBead[®] cell-free DNA Extraction Kit (61C series)

1st Edition, Nov 2022



CONTENTS



Kit description

Introduction	5
Reagent Kits and Instrument	6
Specification	6
Kit Contents	7
Pre-filled Plate or Tube Content	9
Precautions	10
Materials Required but not Provided	10
Kit Storage	11
Safety Information	12
Quality Control	12
Sample Storage	12
Procedure and Description	13

Instrument description

Automatic Nucleic Acid Extraction System	16
Maelstrom 8 (up to 8 tests per run)	17
Maelstrom 4810 (up to 48 tests per run)	18
Maelstrom 9610 (up to 96 tests per run)	19
Instrument Operation Guide	20
Instrument Maintenance	27

Supporting Info

Service Information	28
Technical Support	29
Troubleshooting	30
Ordering Information	33



Kit Description

- 1 Introduction
- 2 Reagent Kits and Instrument
- 3 Specification
- 4 Kit Contents
- 5 Pre-filled Plate or Tube Content
- 6 Precautions
- 7 Materials Required but not Provided
- 8 Kit Storage
- 9 Safety Information
- 10 Quality Control
- 11 Sample Transportation and Storage
- 12 Procedure and Description

Introduction

The TANBead[®] Nucleic Acid Extraction Kit is a nucleic acid purification kit based on magnetic bead technology by using with corresponding TANBead[®] Nucleic Acid Extractor, which can automatically isolate and purify cell-free DNA (cfDNA) from human serum and plasma. The purified cfDNA can be used with many downstream applications employing PCR-based qualitative, semi-quantitative, quantitative assays and capillary electrophoresis. The kit is intended for use by technicians, physicians, and biologists with well-trained in molecular biological techniques, the techniques of magnetic bead purification and in vitro diagnostic procedures. Any diagnostic results generated by using the sample preparation procedure in conjunction with any downstream diagnostic assay should be interpreted related to other clinical or laboratory findings. The kit is not limited to any specific disorder, condition, or other additional accompanying diagnostics. It is applicable for all population.



Reagent Kits and Instrument

Instrument	REF	Cat NO.		Test
		IVD	RUO	
Maelstrom 8 / Maelstrom 4800 series	M61CA46	301385	301508	96
	M61CS46	301389	301509	96
Maelstrom 9600 series	W61CA46	301377	301510	96
	W61CS66	301386	301511	72

Specification

Feature	Specification		
Downstream application	PCR, qPCR and NGS analysis		
Starting material	500-1,000 μ L plasma / serum		
Processing mode	Automated, magnetic processing		
Throughput (Max tests/ run)	Instrument	Auto Plate	Auto Tube
	Maelstrom 8 series	8	4
	Maelstrom 4800 series	48	24
	Maelstrom 9600 series	96	48
Nucleic acid binding technology	Magnetic beads		
Elution Volume	70-100 μ L		
Typical yield	1-100 ng/ per mL plasma or serum		
Typical DNA size	100-300 bp		

Kit Contents

Maelstrom 8 / Maelstrom 4800 series

M61PA46	96 tests	Pre-filled plate for use
Auto Plate	6 PCS	Pre-filled with reagent buffers
Proteinase K	1.2 mL x 2	Proteinase K
Elution Buffer	1.5 mL x 1	Nuclease-free water
Spin Tip	96 PCS x 1 box	Spin tip assembled box
Protocol	1	Instruction guides for user

M61PS46	96 tests	Pre-filled tube for use
Auto Tube	8 Trays	Pre-filled with reagent buffers
Proteinase K	1.2 mL x 2	Proteinase K
Elution Buffer	1.5 mL x 1	Nuclease-free water
Base	2 PCS	A rack for 8 Auto Tubes
Spin Tip	48 PCS x 2 box	Spin tip assembled box
Protocol	1	Instruction guides for user



Kit Contents

Maelstrom 9600 series

M61PA46	96 tests	Pre-filled plate for use
Auto Plate	7 PCS	Pre-filled with reagent buffers
Proteinase K	1.2 mL x 3	Proteinase K
Elution Buffer	1.5 mL x 1	Nuclease-free water
Spin Tip	96 PCS x 1 box	Spin tip assembled box
Protocol	1	Instruction guides for user

W61CS66	72 tests	Pre-filled tube for use
Auto Tube	7 Trays	Pre-filled with reagent buffers
Proteinase K	1.2 mL x 3	Proteinase K
Elution Buffer	1.5 mL x 1	Nuclease-free water
Base	2 PCS	A rack for 8 Auto Tubes
Spin Tip	48 PCS x 2 box	Spin tip assembled box
Protocol	1	Instruction guides for user

Pre-filled Plate or Tube Contents

M61CA46, M61CS46		
Well	Buffer	Volume (μL)
1/7	Lysis Buffer	800
2/8	Washing Buffer 1	1000
3/9	Washing Buffer 2	1000
4/10	Washing Buffer 2	1000
5/11	Magnetic Beads	500
6/12	Elution Buffer	100

W61PC46, W61CS66		
Position	Plate	Volume (μL)
1	Lysis Buffer	800
2	Lysis Buffer	800
3	Washing Buffer 1	1000
4	Washing Buffer 2	1000
5	Washing Buffer 2	1000
6	Magnetic Beads	500
7	Elution Buffer	100
8	Spin Tip	-

Precautions

- 01** Avoid using expired reagents.

- 02** When the temperature is below 20°C, place the Auto Plates / Auto Tubes in an oven (preheated 42~60°C) 5 to 10 minutes.

- 03** Avoid vigorous shaking, to avoid excessive formation of foam.

- 04** Carefully remove aluminum foil to avoid splashing.

- 05** Do not expose the opened reagents or Auto Plates / Auto Tubes to air. The evaporation may cause pH change and affect the extraction effectiveness.

- 06** Please check the integrity of the Auto Plates / Auto Tubes and remember to insert the strips or spin tip into the appropriate position of the suitable instrument before operating them.

- 07** Please wear a mask and disposable gloves when handling.

- 08** Use sterile consumables to avoid nuclease contamination.

- 09** Reagent solution contains guanidine salt, avoid using bleach containing detergent.

- 10** Avoid eyes, skin, and clothing contact with reagents. In case of any contact, flush with flowing water.

- 11** If any serious incident occurs, please report to the manufacturer and the competent authority of the member state in which the user and / or the patient is established.

Materials Required but not Provided

- 01** TANBead® Nucleic Acid Extraction System (non-sterile)
Maelstrom 8/ Maelstrom 4800 series , Maelstrom 9600 series

- 02** Disposable gloves

- 03** Scissors, utility knives

- 04** Micropipette, disposable tips (10 µL / 200 µL / 1000 µL)

- 05** 1.5 mL microcentrifuge tube

- 06** 15 mL / 50 mL conical tube

- 07** 20% SDS (Sodium dodecyl sulfate)

- 08** Streck Cell-Free DNA BCT tubes or compatible ones

Kit Storage

- 01** TANBead® cell-free DNA extraction kit could be stored at 15-35°C for 18 months and please find the expiration date labeled on the box of each reagent kit.

- 02** The proteinase K would be transported at room temperature. Upon received the reagent kit, please keep the proteinase K at 2-8°C for long-term storage. The room temperature delivery of proteinase K has been validated that does not affect its performance.

- 03** The pre-filled plate, tube, and solution provided with the reagent kit can be stored at room temperature and are stable for at least 18 months after production.

Safety Information

When you are working with chemicals, must wear a suitable lab coat, disposable gloves, and safety goggles. For more information, please consult TANBead Inc. for the appropriate safety data sheets (SDSs).

Quality Control

In accordance with TANBead Inc. certificated Quality Management System (QMS and ISO13485), the product must be tested against intended specifications to ensure the consistence of product performance and quality.

Sample Transportation and Storage

01 For short-term storage (up to 7 days), plasma and serum sample should be stored at 2-8°C. For long-term storage, please store samples at -20°C or even -80°C.

02 The storage and transportation condition, such as temperature, time and container may cause the variation of yield and quality of purified nucleic acids. In general, fresh prepared serum and plasma samples are recommended.

03 The transportation of blood samples should follow transportation-related laws and were kept at 2-25°C. Blood samples should be separated for serum or plasma within 6 hrs. The plasma or serum should be transported at 2-8°C or by frozen.

Procedure and Description



TANBead® cell-free DNA extraction procedure containing 3 steps:

01 Sample preparation:

(1) Plasma separation:

- a. The blood sample collected in Streck Cell-Free DNA BCT tubes is stable at 6-37 °C up to 14 days.
- b. Centrifuge at 1,600 x g for 10 minutes.
- c. Transfer the upper layer (plasma) to a new conical tube.
- d. Centrifuge at 16,000 x g for 10 minutes.
- e. Transfer the supernatant (sample) to a new conical tube for cfDNA purification.

(2) Prepare the serum sample:

- a. The blood sample collected in anticoagulant-free collection tubes or serum separator ones.
- b. Keep at room temperature for at appropriate 30 minutes to let clot form.
- c. Centrifuge at 1,600 x g for 10 minutes.
- d. Transfer the upper layer (serum) to a new conical tube.
- e. Centrifuge at 16,000 x g for 10 minutes.
- f. Transfer the supernatant (sample) to a new conical tube for cfDNA purification

02 Sample lysis: Proteins, DNase and RNase can be lysed (inactive) by treating with 20% SDS and proteinase K.



03 Automation: transfer sample to Auto Plate/Tube to proceed lysis, residues washing and elution step for DNA purification.

Plasma sample

STEP 1


Separate the plasma from the blood sample

TIPS
Streck Cell-Free DNA BCT tubes or other cell-free DNA collection tubes are recommended.



STEP 2


After centrifuge, become stratification



STEP 3


Transfer the upper plasma layer to a new conical tube

TIPS
Do not draw the buffy coat as possible to avoid the presence of gDNA in elutes




STEP 4

Carefully transfer the supernatant to a new conical tube for cfDNA isolation



STEP 5

Add 500 μ L supernatant with 30 μ L 20% SDS and 20 μ L Proteinase K (PK)




TIPS
Large sample input (9600 series only)

Sample Volume(μ L)	PK (μ L)	20% SDS(μ L)	Plate Position
500	20	30	Plate 1
1000	40	60	Divide the mixture evenly and add to Plate 1 and 2




STEP 6

Incubate at 60 $^{\circ}$ C for 20 min



STEP 7

Add the total mixture into the well





Serum sample



STEP 1


Keep blood sample at room temperature for ~30 minutes till the clot is appeared.

TIPS
Plasma: Anticoagulant-free collection tubes or serum separator tubes are recommended.



STEP 2



After centrifuge, become stratification



STEP 3


Transfer the upper serum layer to a new conical tube

TIPS
Avoid to draw the clot as possible




STEP 4

Carefully transfer the supernatant to a new conical tube for cfDNA isolation



STEP 5

Add 500 μ L supernatant with 30 μ L 20% SDS and 20 μ L Proteinase K (PK)




TIPS
Large sample input (9600 series only)

Sample Volume(μ L)	PK (μ L)	20% SDS(μ L)	Plate Position
500	20	30	Plate 1
1000	40	60	Divide the mixture evenly and add to Plate 1 and 2




STEP 6

Incubate at 60 $^{\circ}$ C for 20 min



STEP 7

Add the total mixture into the well



Instrument Description



1 Automatic Nucleic Acid Extraction System

Maelstrom 8 (up to 8 tests per run)

Maelstrom 4810 (up to 48 tests per run)

Maelstrom 9610 (up to 96 tests per run)

2 Instrument Operation Guide

3 Instrument Maintenance



Automatic Nucleic Acid Extraction System

Maelstrom 8 (up to 8 tests per run)

Maelstrom 8 handler is a magnetic bead handling device with high-speed stirring function, which can spin up to 3000 rpm. With eight magnetic rods, intuitive interface, and simple operation, it can accomplish any nucleic acid extraction application. It works alone or with Autostage to form an automated solution, called Maelstrom 8 Autostage.



Feature	Specification
REF	Maelstrom 8 handler
Weight (NW)	600 g
Dimensions	11.2(W)x6.3(L)x32.7(H) cm
Power rating	5 Vac, 2A
Battery	3.7 Vac, 2,850 mAh
Max. Throughput	8 samples per run
Process. volume	50 - 1,500 μ L
Spin speed	up to 3,000 rpm
Magnetic rod	> 3,000 gauss
Display	2.4" LCD, 240 x 320 pixels



Maelstrom 4810 (up to 48 tests per run)

Maelstrom 4810 is a 48-throughput instrument, combined with our patented technology, the entire run can be completed in about 15-60 minutes, depending on the reagent kit. Maelstrom 4810 can operate 1 to 48 samples, which offers great flexibility to customers.

Feature	Specification
REF	Maelstrom 4810
Weight (NW)	Approx. 45 kg
Dimensions	58(W) x43(L) x47(H) cm
Power rating	AC 100-AC 240 V 50/60 Hz, 5-2.5 A
Fuse	250 V, 5A
Max. Throughput	48 samples per run
Process. volume	50 - 1,600 µL
Spin speed	up to 3,000 rpm
Heater	12 independent heating blocks
Magnetic rod	> 3,900 gauss
Display	7-inch touchscreen
UV	UV-C type, 8 W
HEPA	E 10 class



Maelstrom 9610 (up to 96 tests per run)

Maelstrom 9610 is a 96-throughput instrument, With the patented technology which can improve the mixing efficiency of magnetic beads and increase the processing sample volume, M9610 has become one of the most competitive automated DNA/RNA extraction instruments. Combined with TANBead extraction reagents, our system is highly affirmed by many medical institutions as it can contribute to large scale sample screening request.

Feature	Specification
REF	Maelstrom 9610
Weight (NW)	Approx. 95 kg
Dimensions	87(W)x57.5(L)70(H) cm
Power rating	AC 220-240 V, 50/60 Hz, 3.5 A AC 100-120 V, 50/60 Hz, 8 A
Fuse	250 V, 5 A
Max. Throughput	96 samples per run
Process. volume	50 - 1,600 µL
Spin speed	up to 3,000 rpm
Heater	4 independent heating blocks
Magnetic rod	> 3,900 gauss
Display	7-inch touchscreen
UV	UV-C type, 4 W
HEPA	E 10 class



Instrument Operation Guide



Maelstrom 4800 series
(Herein represented by Maelstrom 4810)



1


Press the power button in the back of Maelstrom 4800 series instrument



2


2-1


Assemble the spin tip assemble box and select the “Tip” on the panel.




2-2

Choose the intended position, such as 1/7, 2/8, 3/9, 4/10, 5/11, 6/12 and press “Mount”.






Mounting tips (Auto Plate)



Mounting tips (Auto Tube)

3

Carefully remove the aluminum foil on the Auto Plates/Tubes.




4

Add samples to wells of Auto Plate/ Tube, then put it into the extractor.


5

Back to the homepage and choose "Run"



6

Select the corresponding program (program name: 61C)



7


Press "Run" to start the automatic extraction process





Instrument Operation Guide

Maelstrom 9600 series
(Herein represented by Maelstrom 9610)



1


Press the power button in the back of Maelstrom 9600 series instrument




2

Press “Run” and choose the corresponding program
(program name: 61C)





Mounting tips (Auto Plate)




Mounting tips (Auto Tube)

Instrument Maintenance

It is important to clean the device after every use. If samples or reagents have been spilled, it is important to clean the device immediately to avoid damage or contamination of samples.

- 1 Wear gloves and appropriate personal protective equipment.
- 2 If the device is used with biohazardous materials, dispose of any cleaning materials used in accordance with your institutional guidelines.
- 3 The device may go through a run with the magnetic rods unprotected. If this happens, the magnetic rod needs to be cleaned immediately.
- 4 To clean the magnetic rods, wipe with a soft cloth dampened with pure water. Do not use alcohol solvent.

3




Carefully remove the aluminum foil on the Auto Plates/Tubes.

3-1
Add samples to wells of Auto Plate/ Tube.

3-2
Follow the guide on the screen and assemble the required Auto Plate/ Auto Tubes and the assembled spin tip assembled box (NO.8 plate position) to intended plate position.

After placing all plates at right positions, press “Run” to start the automatic extraction process





Supporting Info



1 Service Information

Technical Support

Troubleshooting

2 Ordering Information

Service Information

Technical Support

- 1** Taiwan Advance Nanotech Inc. provides post sales service and technical support in case of any questions, please try to contact our authorized distributor nearest to you or our company for more help.
- 2** If there is a need to contact us, you can call our rep number at +886-3-3167568 or email us at service@tanbead.com for further assistance.
- 3** To efficiently solve the problems and answer your questions, please provide us at least the serial number (for instruments), REF and lot number (for reagent kits) when you talk to our technician personnel.

Troubleshooting

• Less or no nucleic acid yield in the elute

Insufficient sample lysis

Please check the lysis condition (temperature, reagent concentration...etc.) in step ⑥ shown as page 14 and 15. Also, please check the expiry date and storage condition of proteinase K.

Improper pretreatment

Please operate the reagent kit according to the workflow shown as page 14 and 15.

Improper program

Please make sure the extraction program shown on the instrument panel is the same with that shown in IFU.

Sample type and status

The long storage time of sample would have negative effect on its nucleic acid yield, length, quality and integrity. Please find the complete recommendation of storage and transportation shown in section “Sample Transportation and Storage” at page 12.

Troubleshooting



• Magnetic beads remain on the column wall

Improper delivery

This phenomenon is usually caused by transportation and the extraction performance is unaffected. If required, you could centrifuge the plate/tube at 3000 rpm for 1 min or wash them by using the pipetman to make magnetic beads return to the storage buffer.

• Less magnetic beads remain in the elute

Elution contamination

Once magnetic beads remain in all elutes, it may be due to the contaminants in elution buffers, such as microorganism. Please contact our field application engineer (FAE) by phone (+886-3-3167568) or email (service@tanbead.com) for further assistance immediately.

Improper sample input

If this issue appeared in certain samples, please initiate a new extraction again and check whether the issue is still occurred. In most cases of few magnetic beads found in elutes, you could still use the elute for downstream molecular analyses with unaffected extraction performance.

• Amplification failure or unexpected results are found by PCR analysis

Genomic DNA contamination

Please collect separated plasma or serum from blood sample carefully after centrifugation, avoid taking the buffy coat or the clot.

Long amplicon size

Since cfDNA typically shows a length of about 166 bp, for PCR analysis, keep amplicons less than 250 bp is recommended.

Cross-contamination between samples

This issue may be caused by sample overloading. Make sure that the amount of sample and elution buffer keeps within the processing volume.



Maelstrom 8 / Maelstrom 4800 series

Sample	Description	Test	REF	Cat No.
Blood	Blood DNA Auto Plate	96	M611A46	301126
	Blood DNA Auto Tube	96	M611S46	301127
	OptiPure Blood DNA Auto Plate	96	M61EA46	301128
	OptiPure Blood DNA Bulk Plate	960	M61EA10	301307
	OptiPure Blood DNA Auto Tube	96	M61ES46	301129
	Blood RNA Auto Plate	96	M621A46	301400
	Blood RNA Auto Tube	96	M621S46	301401
Plant	Plant DNA Auto Plate	96	M613A46*	301134
		96	M613A46-SE*	301371
	Plant DNA Auto Tube	96	M613S46*	301135
		96	M613S46-SE*	301372
	Plant RNA Auto Plate	96	M6K3A46*	301383
	Plant RNA Auto Tube	96	M6K3S46*	301384
cfDNA	OptiPure cfDNA Auto Plate	96	M61CA46	301385
	OptiPure cfDNA Auto Tube	96	M61CS46	301389
FFPE	OptiPure FFPE DNA Auto Plate	96	M61PA46	301152
	OptiPure FFPE DNA Auto Tube	96	M61PS46	301153
Virus	OptiPure Viral Auto Plate	96	M665A46	301148
	OptiPure Viral Auto Tube	96	M665S46	301149
	OptiPure Viral Bulk Plate	960	M665A10	301346
	Virapid Virus Auto Plate	96	M685A46	301572
	Virapid Virus Auto Tube	96	M685S46	301573
HPV	HPV Auto Plate	96	M61HA46	301589
	HPV Auto Tube	96	M61HS46	301590

Note: black font = IVD and RUO available, blue font = RUO available

Ordering Information

Maelstrom 8 / Maelstrom 4800 series

Sample	Description	Test	REF	Cat No.
Tissue	Tissue DNA Auto Plate	96	M612A46	301130
	Tissue DNA Auto Tube	96	M612S46	301131
	Tissue Total DNA Auto Plate	96	M6T2A46	301132
	Tissue Total DNA Bulk Plate	960	M6T2A10	301306
	Tissue Total DNA Auto Tube	96	M6T2S46	301133
	Tissue Total DNA Auto Kit	96	M6T2046	301260
	Tissue RNA Auto Plate	96	M6K2A46	301366
	Tissue RNA Auto Tube	96	M6K2S46	301367
	Fungi	Fungi DNA Auto Plate	96	M61FA46
Fungi DNA Auto Tube		96	M61FS46	301586
Forensic	Forensic DNA Auto Plate	96	M6TFA46	301424
	Forensic DNA Auto Tube	96	M6TFS46	301425
Bacteria	Gram Bacteria DNA Auto Kit	96	M61G046	301257
	Gram Bacteria DNA Auto Plate	96	M61GA46	301138
		96	M61GA46-SE	301294
	Gram Bacteria DNA Auto Tube	96	M61GS46	301139
		96	M61GS46-SE	301295
Plasmid	Plasmid Extraction Auto Plate	96	M6PEA46*	301578
	Plasmid Extraction Auto Tube	96	M6PES46*	301579
Food Feed	Food and Feed DNA Auto Plate	96	M6GMA46*	301635
	Food and Feed DNA Auto Tube	96	M6GMS46*	301636
Environmental Microbiome	Environmental Microbiome DNA Auto Plate	96	M6EMA46*	301641
	Environmental Microbiome DNA Auto Tube	96	M6EMS46*	301642
Stool	Stool Cell DNA Auto Plate	96	M6SCA46	301387
	Stool Cell DNA Auto Tube	96	M6SCS46	301388

Note: black font = IVD and RUO available, blue font = RUO available

Ordering Information



Maelstrom 9600 series

Sample	Description	Test	REF	Cat No.
Blood	Blood DNA Auto Plate	96	W611A46	301186
	Blood DNA Auto Tube	72	W611S66	301187
	OptiPure Blood DNA Auto Plate	96	W61EA46	301188
	OptiPure Blood DNA Auto Tube	72	W61ES66	301189
	Blood RNA Auto Plate	96	W621A46	301402
	Blood RNA Auto Tube	72	W621S66	301403
	Plant	Plant DNA Auto Plate	96	W613A46*
96			W613A46-SE*	301379
Plant DNA Auto Tube		72	W613S66*	301259
		72	W613S66-SE*	301378
Plant RNA Auto Plate		96	W6K3A46*	301406
Plant RNA Auto Tube		72	W6K3S66*	301407
cfDNA		OptiPure cfDNA Auto Plate	96	W61CA46
	OptiPure cfDNA Auto Tube	72	W61CS66	301386
FFPE	OptiPure FFPE DNA Auto Plate	96	W61PA46	301629
	OptiPure FFPE DNA Auto Tube	72	W61PS66	301630
Virus	OptiPure Viral Auto Plate	96	W665A46	301224
	OptiPure Viral Bulk Plate	960	W665A10	301345
	OptiPure Viral Auto Tube	72	W665S66	301209
	Virapid Virus Auto Plate	96	W685A46	301574
	Virapid Virus Auto Tube	72	W685S66	301575
HPV	HPV DNA Auto Plate	96	W61HA46	301591
	HPV DNA Auto Tube	72	W61HS66	301592

Note: black font = IVD and RUO available, blue font = RUO available

Ordering Information








Maelstrom 9600 series



Sample	Description	Test	REF	Cat No.
Tissue	Tissue DNA Auto Plate	96	W612A46	301190
	Tissue DNA Auto Tube	72	W612S66	301191
	Tissue Total DNA Auto Plate	96	W6T2A46	301192
	Tissue Total DNA Auto Tube	72	W6T2S66	301193
	Tissue RNA Auto Plate	96	W6K2A46	301404
	Tissue RNA Auto Tube	72	W6K2S66	301405
Forensic	Forensic DNA Auto Plate	96	W6TFA46	301291
	Forensic DNA Auto Tube	72	W6TFS66	301426
Bacteria	Gram Bacteria DNA Auto Plate	96	W61GA46	301198
	Gram Bacteria DNA Auto Tube	72	W61GS66	301199
Plasmid	Plasmid Extraction Auto Plate	96	W6PEA46*	301580
	Plasmid Extraction Auto Tube	72	W6PES66*	301581
Food Feed	Food and Feed DNA Auto Plate	96	W6GMA46*	301637
	Food and Feed DNA Auto Tube	72	W6GMS66*	301638
Environmental Microbiome	Environmental Microbiome DNA Auto Plate	96	W6EMA46*	301643
	Environmental Microbiome DNA Auto Tube	72	W6EMS66*	301644
Stool	Stool Cell DNA Auto Plate	96	W6SCA46	301392
	Stool Cell DNA Auto Tube	72	W6SCS66	301391

Note: black font = IVD and RUO available, blue font = RUO available

Consumables




Sample	Description	Test	REF	
	96 deep well plate	Auto Plate	<ul style="list-style-type: none"> 100 pcs/carton Processing volume 50-1,600 µL For molecular diagnostics 	083.MWP01.20X
	96 deep well plate (Unique hook design)	Auto Plate	<ul style="list-style-type: none"> 100 pcs/carton Processing volume 50-1,600 µL For molecular diagnostics 	083.MWP02.20X
	Spin Tips Assembled Box (Unique hook design)	Auto Plate	<ul style="list-style-type: none"> 80 pcs/carton 96 pcs of spin tips in one box 	083.MSP09.10X
	Spin Tips Assembled Box (Unique hook design)	Auto Tube	<ul style="list-style-type: none"> 80 pcs/carton 48 pcs of spin tips in one box 	083.MSP10.10X
	6 tube B	Auto Tube	<ul style="list-style-type: none"> 96 pcs/bag, 16 bags/ cartonme Minimal consumable waste 	104143
	16-Base B	Auto Tube	<ul style="list-style-type: none"> 300 pcs/carton Integrate with 6 tube B for small number of tests 	104026
	Spin tips	Auto Plate Auto Tube	<ul style="list-style-type: none"> 2,000 pcs/bag, 10 bags /carton A unique design for maximum mixing efficiency 	056.CSM03.111





Taiwan Advanced Nanotech Inc.

 6F., No. 188, Wenhe Rd., Guishan Dist.,
Taoyuan City 333, Taiwan

 TEL / +886-3-3167568
FAX / +886-3-3173369

 success@tanbead.com

 <https://www.tanbead.com>

